*Azura*

**Mixing chamber**

**Mischkammer**

Supplement/Beiblatt

**HPLC**

Document no. V6730

**Note:** For your own safety, read the corresponding instructions and observe the warnings and safety information on the device and in the instructions.

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Version number: 1.3

Release date: May 05, 2023

The information in this document is subject to change without prior notice. For the latest version of the instructions, visit our website: [www.knauer.net/library](http://www.knauer.net/library).

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# 1. Mounting the Mixer

You mount the mixer with 2 Allen screws. The screws are included with the accessory kit.

- Tools:**
- Allen screwdriver, 2.5 mm
  - torque wrench, 1/4"

**Note:** The biocompatible mixer can be used with stainless steel or PEEK fittings for the capillaries. In case you use stainless steel, the sealing ring must be made from PEEK.

## NOTICE


### Component defect

Damage to components due to excessive tightening possible.  
Observe the torque of the screw connection.

- Use 5 Nm torque for stainless steel fittings.
- Use 1 Nm torque for PEEK fittings.

## Process and Figure

1. Attach the mixer and fix it with the screws ①.
2. Using the screwdriver, fasten the screws.
3. Fasten the screw fittings ② together with the capillaries.

**Fig. 1** Mounting the mixer

Then put the device into operation.

## 2. Dismounting the Mixer


The mixer is dismounted in one piece, you are not allowed to disassemble the subassembly.

- Prerequisites:**
- The mixer has been flushed with isopropanol.
  - The hole plugs are laid out ready.

- Tools:**
- open-end wrench, size 1/4"
  - Allen screwdriver, 2.5 mm

### Process

1. Disassemble the screw fittings ② by hand or with the open-end wrench.
2. Close the openings with the hole plugs.
3. Using the screwdriver, unfasten the screws ①.
4. Remove the mixer and lay it to the side.


**Fig. 2** Dismounting the mixer

Connect the capillaries, so there is no way for solvent to leak. Then put the device into operation.

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All rights reserved.
AZURA® is a registered trademark of
Wissenschaftliche Geräte GmbH