Smartline UV Detector 2600

Manual / Handbuch

V 7605 05/2007

Smartline

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Contents

Using this Manual	7
Conventions in this manual.....	7
SOP's in this manual	7
General Description of Detector 2600	8
General Description of Detector 2600	8
The Basic Principle of the Detector 2600	9
Optical Path of the Detector 2600.....	9
Preparing the Detector 2600 for Operation.....	10
Unpacking	10
Standard Delivery Smartline Detector 2600	10
Front and Rear View of the Detector 2600	11
Front Door of the Smartline Detector 2600.....	11
The Metal Front Plate of the Smartline Detector 2600	11
Rear Panel of the Smartline Detector 2600.....	12
Power supply, ON/OFF	13
The Detector's Position in a KNAUER Smartline system	13
Operating the Detector 2600.....	14
Function of Foil Keys	14
AUTOZERO key	14
SCAN key	14
Numeric keys	14
Arrow keys	14
Standby key	14
Switching on the Detector	15
Display and Structure of the Menus.....	15
Main Menu	15
Appearance of the Cursor and Data Input.....	16
SETUP Menu	17
Time constant	17
Bandwidth	18
Lamp	18
Signal	18
Scan out.....	18
Analog out.....	18
Analog out offsets	19
Integration time	19
Date and time	19
Fraction level	19

Fraction event 1 and 2	19
Fraction valve.....	19
Shutter.....	19
LCD Contrast	20
GLP menu	20
Lamp	21
Total working time lamp	21
Total working time device	21
Last error codes	21
Installation date	21
Calibration date	21
Last service	21
Serial number.....	21
Wavelength validation.....	21
Wavelength validation results	22
Firmware Version.....	22
VIEW menu	22
LINK menu	23
Installation of the Flow Cell	23
Capillary Connection to an HPLC System	26
Setting the Integration Time	27
Direct Control of the Detector 2600 (stand-alone mode).....	28
Wavelength Selection	28
EVENTS.....	29
Programming the Smartline Detector 2600.....	30
Creating Programs.....	30
Event Programming	32
Fraction Collector Control	32
Program execution.....	35
Creating Program Links	35
Programming Links	35
Running Program Links	36
Deleting Programs and Program Links.....	36
Wake Up Program	37
Using the Wake Up Program	37
Scanning UV Spectra.....	38
The Scan Menu.....	38
Performing Wavelength Scans	38
Output of Scan Data	38
Connecting other Instruments to the Detector 2600	39
Using the Remote Sockets.....	39

Connecting the Remote Sockets	39
Assembling WAGO Plugs.....	40
The Integrator Outputs.....	41
RS-232 Serial Interface	41
USB Interface	42
Software Control of the Detector 2600	42
Simple Maintenance	43
Control of the Lamp functionality	43
Checking Wavelength Accuracy	43
Changing the lamp.....	45
Cleaning the flow cell.....	47
Analytical flow cells.....	48
Analytical 3 mm Flow Cells.....	48
Preparative Flow Cells.....	49
Adjusting the Path Length of Preparative Flow Cells.....	50
Error Messages.....	51
Error Messages, Descriptions and Codes	51
Error Messages and their Effects	51
Error Messages and their Reasons	52
Other Problems and possible Causes	53
Spare Parts and Accessories.....	54
Flow cells for the Smartline Detector 2600.....	54
Analytical Flow Cells.....	54
Preparative Flow Cells.....	54
U-Z View™ Micro Flow Cells.....	54
Spare Parts and Accessories, Order Numbers	54
Specifications	55
Warranty Statement.....	107
Declaration of conformity	108
Index	109

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Using this Manual

This manual refers to the Smartline Detector 2600 (UV or VIS) Firmware Revision 1.03 or higher. It is valid for any combination with analytical flow cells order number A4061, A4042, A4045, preparative flow cells order number A4066, A4067, A4069 and A4095.

Conventions in this manual

Arrows like these: ◀ ▶, used in block diagrams, indicate that the user is asked to press the corresponding arrow keys. The operation of arrow keys is defined as follows:

Cursor up: ▲

Cursor left: ◀

Cursor right: ▶

Cursor down: ▼

Important Hints are marked by the marginal hand symbol.

Special warnings are indicated by the marginal warning sign and printed in bold letters.

The marginal lamp symbol indicates helpful advice.

SOP's in this manual

The **Standard Operating Procedures (SOP)** provided with this manual offer a convenient way of structuring complex tasks in the operation of your **Smartline Detector 2600**. They include step-by-step instructions leading the user through all routine tasks during operation. They can be used for documentation purposes and be copied, applied signed, and filed in order to document the performance of the instrument.

Please operate the instrument and all accessories according to instructions and SOP's in this manual. This ensures safe operation, accurate results and prolongs the longevity of your equipment.

SOP 1	Installation of the Flow Cell	24
SOP 2	Capillary connections	26
SOP 3	Automatic Integration Time Determination	27
SOP 4	Time Programming of Wavelengths	30
SOP 5	Event Programming for Fraction Collector Control	33
SOP 6	Assembling WAGO plugs	41
SOP 7	Wavelength accuracy check	44
SOP 8	Changing the deuterium lamp.....	45
SOP 9	Changing the tungsten halogen lamp.....	47
SOP 10	Purging the Flow Cell	47
SOP 11	Disassembling and Cleaning of an Analytical 10 mm Flow Cell.....	48
SOP 12	Disassembling and Cleaning of an Analytical 3 mm Flow Cell	49
SOP 13	Disassembling and Cleaning of a Preparative Flow Cell.....	49
SOP 14	Changing the Path Length of Preparative Flow Cells	50

General Description of Detector 2600

Fig. 1 Smartline Detector 2600

The Smartline Detector 2600 is a completely new multiwavelength spectrophotometer that is based on diode array technology. The instrument unites the advantages and the price of an outstanding 4-channel spectrophotometer with the quality of a diode array detection system. The result is a detector with unique performance and flexibility in the new compact KNAUER Smartline design.

The Smartline Detector 2600 can simultaneously measure at up to four freely programmable detection channels with high sensitivity and continuously record spectra in a selectable scan range. The signals of the four channels are additionally distributed on the four analog outputs. The instrument operates depending on the installed lamp in a wavelength range of 190 to 510 nm (UV version) or 430 to 710 nm (VIS version) with a selectable bandwidth of 4 to 25 nm. Calibration and wavelength validation are supported by an integrated holmium filter.

The flow cell is located behind the front panel, on which the display and touchpad are integrated. It is well protected from outside influences, yet made easily accessible. All current KNAUER flow cells can be used in the Smartline Detector 2600.

The entire optical unit can be removed from the front of the detector, making lamp exchanges easy since opening the housing is not necessary. Due to the fact that the instrument has practically no moving parts, it requires very little maintenance and is highly robust.

The Smartline Detector 2600, together with KNAUER ChromGate®, ClarityChrom®, or EuroChrom® Software, enables all functions that are typical of a diode array detector including multi-wavelength operation, 2-D and 3-D chromatograms, peak purity, spectra libraries, and more. In addition, control is possible through the remote control connection, the open communication protocol, or the touchpad (stand-alone operation).

Smartline Detector 2600 – Everything is possible.

The Basic Principle of the Detector 2600

Optical Path of the Detector 2600


Fig. 2 Optical path of the Detector 2600 (diagram)

The light from the deuterium lamp (1) is focused by the collimator lens (2), passes the filter wheel (3) and then goes through the flow cell (4). Two mirrors (5) direct the beam in the direction of the condenser lens (6) which projects the light on the slit (7). After passing the slit, the light is reflected by the first concave mirror (9) onto the grating (8) where the whole spectrum is split up and then reflected by the second concave mirror (9) on the photodiode array (10). This photodiode array is placed in a way that the first diode represents 190 nm or 430 nm and the last one represents 510 nm or 710 nm. Every single photodiode measures the absorption at the corresponding wavelength. The Smartline Detector 2600 needs no reference beam.


Fig. 3 Optical path of the Detector 2600

Preparing the Detector 2600 for Operation

Unpacking


All KNAUER instruments are always packed carefully and safely for transportation. After unpacking, please check the device and accessories thoroughly for any damage that may have occurred during transport. If necessary, put forward any claim for damages to the carrier.

Use the list "Standard delivery" and check if the Detector 2600 delivery is complete. Please contact our service department if you are missing something or if you need support.


Remove the transport protections from the display and from the flow cell housing.

Standard Delivery Smartline Detector 2600


1. Smartline Detector 2600 with test cell, without flow cell


2. Operation manual

3. Power supply cable
4. RS-232 Cable (9 pin, female/female)
5. four analog connection cables (cinch/cinch)
6. two WAGO connectors (8-fold) and mounting device
7. two flat cables (10 conductor)
8. Allen key (size 3)

9. Holmium filter certificate

Front and Rear View of the Detector 2600

Front Door of the Smartline Detector 2600


Fig. 4 Front view of the Detector 2600

1	Display
2	Arrow keys (up, down, left, right)
3	Standby key
4	Auto-Zero key
5	Scan key
6	Numeric keys

Table 1 Functional elements on the front

The Metal Front Plate of the Smartline Detector 2600

(behind the door)


Fig. 5 Front view with open door

1	Knurled screws that hold the flow cell
2	Screw for removing the optical bench (size 3)

Table 2 Functional elements on the front plate


Fig. 6 analytical flow cell

1	Inlet of the cell
2	Outlet of the cell
3	Indication of the path length (10 mm)

Table 3 analytical flow cell

Rear Panel of the Smartline Detector 2600

Fig. 7 Rear view of the Detector 2600

1	Analog outputs (1V, scalable)
2	USB (not active !)
3	Remote control connectors 1
4	RS-232 connector
5	Remote control connectors 2
6	Serial number
7	Power connector
8	Main switch

Table 4 Functional elements on the rear panel

Power supply, ON/OFF

The Detector 2600 is equipped with a modern switching power supply. This power supply works in a range of 90 - 260 V AC and 47-63 Hz. The instrument can be switched off from the main switch (on the back of the instrument) or the standby key.

Please note that the instrument is not completely switched off when in standby-mode. Only with the main switch the instrument can be completely switched off.

The Detector's Position in a KNAUER Smartline system

Due to the general temperature sensitivity of detectors, the Smartline Detector 2600 should always be the first (lowest) instrument in a Smartline system. The optionally available capillary kits for easy installation will only fit if this rule is followed.

Fig. 8 shows a typical Smartline system (low pressure gradient).


Fig. 8 Smartline system (low pressure gradient)

In the case that a Smartline RI Detector should also be arranged in the tower, the RI detector should be installed below the UV or VIS detector due to its higher sensitivity to temperature.

Operating the Detector 2600

Function of Foil Keys

The keypad in Fig. 4 consists of ten numeric keys, an AUTOZERO key, a SCAN key, four arrow keys and a standby key.

AUTOZERO key

Pushing this button will adjust the baseline. Usually the button should be pressed shortly before starting a chromatographic run. Pressing the AUTOZERO key will set all acquisition channels to zero.

SCAN key

This key is used for manual spectra acquisition. After pressing the SCAN key, the SCAN menu is shown. In this menu you can either select between scan nos. 1-4 or you can send a recorded spectrum to analog output 1. You can also record a scan by pressing the SCAN key again. The individual wavelength range of the detector is also shown in this menu.

Numeric keys


The numeric keys are used to enter numerical values at the recent cursor position (if possible).

In certain cases the numeric keys are used e.g. to start (1) or interrupt (0) a program.

Arrow keys

The arrow keys in light green are used for moving the cursor on the display and to confirm settings.


Use the arrow keys \triangleright "right" or \triangleleft "left" to move from one menu parameter to the other as well as to confirm entered or selected values.


Using the arrow keys \triangleup "up" or \triangledown "down" you can change the selected parameter or the available options of it.

Standby key

Pressing the standby key for more than two seconds switches the instrument off (only the standby electronics will still be working). In standby mode the red LED integrated into the standby key is on. To switch the instrument on again the standby key has to be pressed for more than one second. Then the instrument powers up again and the red LED is off.

Switching on the Detector

Connect the detector to the power supply cable and switch on the instrument. The main power switch is on the rear panel.

Detector 2600

After switching ON, the display presents briefly information concerning the instrument and the firmware version:

**Firmware Version 1.01
Build 1 14.01.04**

Whenever powering up, the instrument includes initialization and check of the electronics and a check of the filter wheel. Also the lamp is heated to its working temperature for about 20 seconds and then it is switched on. The whole procedure is indicated by the following display:

SYSTEM TEST

After successfully completing the self test the detector is ready for use. The main menu is shown on the display. Heating of the lamp is indicated by 'HEA' and the message 'ON' shows that the lamp was switched on.

```
-- -- 000 .3219 .---
ON 000.00 254 000
```

After a period of about 10 to 30 minutes the instrument reaches a constant working temperature. The baseline is now stabilized and the UV or VIS detector respectively is ready to acquire data.

Display and Structure of the Menus

After switching on the detector the instrument is ready for operation and the main menu is shown on the display.

The blinking cursor can be moved to all positions in the menu by using the arrow keys. If the cursor reaches a position at the edge of a line (e.g. **ON** in the lower left corner of the main menu) then it is possible to enter the corresponding menu (here: the VIEW menu) by pressing the appropriate arrow key (here: the \triangleleft key) again for about two seconds. There are four additional menus: The LINK menu, the VIEW menu, the SETUP menu and the GLP menu. To return to the main menu you can press the reverse arrow key (here: the \triangleright key) one or several times depending on the position in the additional menu. If no key was pressed for more than 60 seconds the instrument automatically returns to the main menu.

Main Menu

The main menu displays the signal values (**4/5**) and wavelengths (**8/9**) for channel 1 and channel 2 (λ_1/λ_2). Furthermore there is a field for loading programs (**1**), a field for editing programs (**2**), two fields that show the current state of the event outputs (**3**) and the current state of the lamp (**6**) and a field that shows the current time (**7**).

SETUP Menu


Fig. 12 SETUP menu overview

Parameters for the operation of the Smartline Detector 2600 can be adjusted from the SETUP menu. These parameters comprise all of the necessary values for operation of the instrument.

Time constant

TIME CONSTANT: 1.0s

The signal can be smoothed with a time filter. The time constant value can be chosen between **0.1**, **0.2**, **0.5**, **1**, **2**, **5** and **10** seconds. The larger the value for the time constant, the more the signal is smoothed. For most analytical purposes, a time constant of **1s** is most suitable.

Bandwidth

BANDWIDTH: 04 nm
 ◆ Range: 04-25 nm

The bandwidth can be selected between 4 and 25 nm. Different bandwidths can be implemented by 'diode bunching'. The lower limit (4 nm) is given by the real optical bandwidth of the detector.

Lamp

LAMP: Deuterium ON
 ◆ no: 002 00095.7

The lamp type is set to **deuterium** and cannot be changed. The current state of the lamp (ON/HEA/OFF) is shown together with the number of lamps used in this instrument and the operation time of the current lamp (hours).

Here you can switch the lamp on and off (**ON/OFF** field) and in case the lamp was replaced the lamp counter can be increased (no:). At the same time, the operation time will be reset.

Please increase the lamp counter only after replacing the lamp ! The lamp counter can only be increased. It can not be edited. In addition, the lamp counter will not reset from '999' (maximum value) back to '000'.

Signal

SIGNAL: sig1:0.4544
 ◆ absorption sig2:0.0000

SIGNAL: sig3:0.0000
 ◆ absorption sig4:0.0000

The signal can be shown as either absorption or intensity (absorption / sig-channel). All channels together (max. four) are shown as absorption or intensity. It is not possible to select the mode for each channel separately.

Scan out

SCAN OUT
 ◆ speed: 1 nm/s

Stored spectra can be put out via analog output 1 (analog outputs 2, 3 and 4 are not used for spectra output). The speed for spectra output can be set to either 1, 2, 5, 10, 20 or 50 nm/s.

Analog out

ANALOG OUT FULL SCALE:
 ◆ 1 au-> 1V

The maximum signal of the four analog outputs is 2,5 V. In this submenu you can choose which absorption level should equal 1 V. You can choose 0,0001, 0,0002, 0,0005, 0,001, 0,002, 0,005, 0,01, 0,02, 0,05, 0,1, 0,2, 0,5 1, 2, 5 or 10 AU. Additionally the voltage can be changed from 1 V to 0.1 V.

Analog out offsets

ANALOG OUT OFFSETS 1 2:
 ◆ +000.00 mV +000.00 mV

ANALOG OUT OFFSETS 3 4:
 ◆ +000.00 mV +000.00 mV

For each analog output an individual offset (between -999.99 mV and 999.99 mV) can be defined. You can switch between + and – with the arrow keys ‘up’ (▲) and ‘down’ (▼).

Integration time

INTEGRATION TIME: 020 ms
 ◆ Range: 13-200 ms

The integration time can be set manually or automatically between 13 and 200 ms. Default value is 20 ms. In standalone operation the optimum integration time should be set automatically before data acquisition. Under software control setting the integration time is not necessary because this is done by the software.

Date and time

DATE: 12.02.04 dd.mm.yy
 ◆ **TIME: 16.39.30 hh.mm.ss**

Date and time will be shown as daymonth, monthmonth, yearyear and hourhour, minuteminute, secondsecond.

Fraction level

FRACTION LEVEL
 ◆ 0.00000 AU

If the fraction collection feature is to be used, a certain minimum absorption level must be defined. At or above this level, the fraction events (or fraction valve) will be activated. This level can be selected between 0.00001 AU (10 µAU) and 9.99999 AU.

Fraction event 1 and 2

FRACTION EVENT1
 ◆ OFF DELAY: 000.00 sec

FRACTION EVENT2
 ◆ OFF DELAY: 000.00 sec

Simple fraction collection can be realised using the optional fraction events 1 and 2 (which affect event outputs 1 and 2) as well as the fraction valve output. A detailed description of the fraction collection functionality can be found on page 32.

Each fraction event (or fraction valve) needs to be activated or deactivated (ON ↔ OFF) separately. If necessary, an individual delay (0.01-999.99 s) can be added.

Fraction valve

FRACTION VALVE
 ◆ OFF DELAY: 000.00 sec

The fraction collection valve can be activated by switching from OFF to ON. If necessary an individual delay (0.01-999.99 s) can be added.

Shutter

SHUTTER: OPENED
 ◆

Shutter refers to the filter wheel. There are three positions available (CLOSED/OPENED/FILTER). CLOSED means that no light gets to the cell and subsequently to the diode array. In this position the dark current can be measured. OPEN is the standard setting and is used for normal data acquisition. FILTER means that a holmium oxide filter glass is moved into the optical path. This filter is used to determine wavelength accuracy.

LCD Contrast

LCD CONTRAST: 00
◆

The Display's contrast can be set between 00 to 10 (00 = dark, 10 = light).

GLP menu

The GLP menu provides detailed operation and service data. GLP data can not be edited. In addition the automatic wavelength validation can be run here.


Fig. 13 GLP menu overview (in case of a tungsten halgen lamp the HO λ values 637,7; 536,5; 445,7 and 400,2 will be displayed)

Lamp

GLP: Lamp no:0001
◆ 067.2 h starts:0018

This submenu shows the number of lamps used in this instrument (no:), the operation time of the current lamp and the current number of lamp starts (starts:).

Total working time lamp

GLP: Total working time
◆ (lamp) 0067.2 h

Displays the total working time of the current lamp.

Total working time device

GLP: Total working time
◆ (device) 0067.2 h

Displays the total working time of the instrument.

Last error codes

GLP: Last error codes
◆ 0 0 0 0

Shows the error codes of the last four errors. The Table of error codes can be found on page 51.

Installation date

GLP: Installation date
◆ 05.02.04

Shows the day on which the instrument was built.

Calibration date

GLP: Calibration date:
◆ 11.02.04

Displays the date of the last calibration.

Last service

GLP: Last service
◆ 00.00.00 code:0000

Displays the last service date.

Serial number

GLP: Serial number:
◆ 12345

Shows the instrument's serial number.

Wavelength validation

Wavelength validation
◆

From this menu the automatic wavelength validation can be accessed. The wavelength validation is running fully automatically using the

detector's integrated holmium oxide glass filter. Altogether four lines of the holmium filter will be checked. If the deviations from the stored line positions are ≤ 1 nm the validation is passed. In case one or more deviations is/are higher the validation is failed. All results are stored until the next validation is performed. If the instrument failed in the validation the according error code (error code 3) is shown in the list of last errors.

Wavelength validation results

HO λ	Meas. λ	1
◆ 241.5	240.8	OK

HO λ	Meas. λ	2
◆ 279.2	278.8	OK

HO λ	Meas. λ	3
◆ 361.0	360.6	OK

HO λ	Meas. λ	4
◆ 446.0	445.4	OK

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HO λ	Meas. λ	1
◆ 240.9	240.6	OK

HO λ	Meas. λ	2
◆ 445.7	445.3	OK

HO λ	Meas. λ	3
◆ 536.5	536.1	OK

HO λ	Meas. λ	4
◆ 637.7	637.2	OK

Tungsten halogen lamp

The results of the last wavelength validation are shown. For four holmium lines the real wavelength (HO λ), the measured wavelength (Meas. λ) and the result (OK or Failed) are displayed.

Firmware Version

GLP: Firmware Ver. 1.01
◆ Build 1 14.01.04

Displays the installed firmware version.

VIEW menu

The View menu gives an overview of existing programs and program links. See Fig. 14 for an example where programs 1, 2, 4 and 12 and link 26 have been defined. Program no. 30 is the WAKE UP program.

Main Menu

--- 000 .4517 ---
◆ N 000.00 λ 3 300 λ 4 000

VIEW Menu

Programs 0 1 02 --- 04 ---
◆ (01..10) --- --- --- ---

△ ▽

Programs:--- 12 --- ---
◆ (11..20) --- --- --- ---

Programs:△ ▽ --- --- --- ---
◆ (21..29) 26 --- --- --- 30

Programs 0 1 02 --- 04 ---
◆ (01..10) --- --- --- ---

Fig. 14 VIEW menu overview

LINK menu

The link menu is used when running links, i.e. the predefined run of program combinations. This menu shows the current state of the links (Fig. 16). If no link is active only the standard display is shown (Fig. 15).


Fig. 15 LINK menu without active link

shows the LINK menu during operation. Link L1 (program no. 21) is running. The currently running program (No01) is shown in the first line. The program which is to run next is shown in the second line. P01 refers the number of the currently running program and R-05 means five repetitions.

First only two program lines (No01 and No02) can be seen. P01 means execution of program 01 and R05 means 5 repetitions (of program No 01). The second line says that afterwards program no. 13 will be repeated 10 times. Using the ∇ key it is possible to move further down to line no. 2 (No02) and 3 (No03).

For more information about links, see the chapter about programming and running links on page 35.


Fig. 16 LINK menu with active link (program line No01 is running)

Installation of the Flow Cell

The Detector 2600 comes factory configured without a flow cell. The device is equipped with a 'dummy' test cell which does not have any optical parts and is for testing the instrument only. Before operating the detector it is necessary to install an appropriate KNAUER flow cell.

A list of available flow cells can be found on page 54.


Fig. 17 Test cell ("dummy" cell), left; analytical flow cell, right

SOP 1 Installation of the Flow Cell

1. Swing open the detector's front door.
2. Loosen and remove the two knurled screws (1) on the cell holder's front (2).
3. Pull out the cell holder.
4. Hold the cell or dummy cell (3) with two fingers and remove it upwards.
5. Insert the new flow cell (5) and make sure that the outlet is on the upper side and that the fixing hole on the back side of the cell meets the corresponding metal pin (4) of the cell holder.
6. Now push the complete system towards the housing, insert the two screws (1) and tighten them manually.

The instrument and lamp need not be switched off when exchanging the flow cell. The instrument is ready for immediate operation after proper installation of a flow cell.


Fig. 18 Deinstalling the test cell (dummy cell)


Fig. 19 Deinstalling the test cell (3) and installing a flow cell (5)

Capillary Connection to an HPLC System

Before taking a measurement cell filled with solvent into operation, please make certain that the eluent used is miscible with the one used previously. Otherwise, purge the flow cell with a medium which is miscible with both of the eluents.

Even though the Smartline Detector 2600 is very resistant to many kinds of commonly used eluents, you should take care that no eluent or water can get on the instrument's surface or even inside the instrument. Chlorinated hydrocarbons could destroy the varnish and some others (e.g. THF) could loosen the keypad, for example.

SOP 2 Capillary connections

1. Connect the outlet of the HPLC column to the inlet bushing of the flow cell (lower side).
2. Push the bushing, the clamping ring, and the sealing ring onto the capillary (for standard connections just use the respective bushing and clamping ring). Please take care to note the sequence and orientation of the fittings according to Fig. 21
3. Push the capillary as far as possible into the flow cell inlet.
4. Fasten the bushing by hand (for standard stainless steel connection: please use the supplied wrench).
5. Connect the flow cell outlet (upper side) using a capillary or PTFE tube (ID > 0.5 mm) to a waste bottle.

The capillary connections in a simple HPLC system are shown in the following figure.

Use bushings (e.g. DYNASEAL bushings, see Accessories on page 54) which keep the dead volume as small as possible and the shortest possible capillary with a small internal diameter.


Fig. 20 Capillary connections of the detector
(example: Low pressure gradient HPLC system)


Fig. 21 DYNASEAL Capillary connections

Setting the Integration Time

Setting the integration time is necessary in standalone mode only.

Automatic determination of the integration time is useful only if it is performed with the recent flow cell and eluent.

SOP 3 Automatic Integration Time Determination

1. Move the cursor to the INTEGRATION TIME sub menu in the SETUP menu.

INTEGRATION TIME: 020 ms
◆ Range: 13-200 ms

- Press the \triangleright key and enter a value between 13 and 200 in the **INTEGRATION TIME** field. Alternatively you can press the \triangle key and the optimum integration time will be determined automatically.

INTEGRATION TIME: 20 ms
 \blacklozenge cursor up = aut.define

- During the automatic determination the following message is shown:

INTEGRATION TIME: 020 ms
 \blacklozenge please wait...

If the integration time was set too low the signal/noise ratio would be low. On the other hand a very high integration time may lead to saturation in areas with high light intensity (e.g. around 240 nm). In areas which show saturation no signal will be detected.

Direct Control of the Detector 2600 (stand-alone mode)

Switch the instrument on, keeping in mind the hints given in the section Power supply, ON/OFF on page 13.

Prior to the first measurement wait for about 15 minutes for the instrument to warm up while the HPLC pump is also switched on (with flow). In case of especially sensitive measurements prolonging this warming up period may be necessary.

The lamp starts automatically whenever the detector is switched on. However, at any time it can be switched off and on in the SETUP menu (Lamp sub menu) (see page 18).

Set the desired wavelength in the main menu and then change to the SETUP menu. There you can select the appropriate time constant and other settings.

In principle, now your Detector 2600 is ready to take simple chromatograms.

The Smartline Detector 2600 can be controlled directly via the front keypad. Numerous functions like for example setting the wavelength (except scanning) can be chosen by using numeric and arrow keys.

Programs can be used for wavelength control and for setting events during a run.

This Smartline Detector 2600 offers up to 20 programs with up to 100 lines (total). In addition the instrument offers 9 link programs with up to 100 lines. Link programs are used to link different programs together in one large link program.

Due to the working principle of the detector (diode array) it can record up to 10 spectra per second. Scans are executed from the SCAN menu which is easily accessible via the SCAN key.

Wavelength Selection

Position the cursor on one of the two wavelength fields and select the wavelength by entering the desired values with the numeric keys. Wavelength three and four (λ_3 und λ_4) can be entered after moving the cursor further to the right. If you want to use only one single wavelength,

enter just nothing or a zero in the unused wavelength fields. Entering zero deactivates the respective channel. If the detector works at two or more wavelengths this has no influence on the detector's noise as the instrument performs scans the whole range from 190 to 510 nm (UV version) or 430 to 710 nm (VIS version) anyway.

Output options: The output of the absorption signals can be configured. Move the cursor to the signal fields or to the wavelength fields. Choose between different options by pressing any numeric key. Options are:

- 1) **Normal signal** (normal presentation)
- 2) **Signal negation** (indicated by a little superimposed minus sign). Signal negation for the corresponding channel is activated and deactivated by pressing any numeric key when the cursor is located on a signal field.
- 3) **Ratio plot 1** (the symbol λ_2^1 is shown above the decimal point of the absorption signal) The display shows the result of the division of signal (channel 1) by signal (channel 2).
Ratio plots can be selected on channel 3 and 4. They are activated, changed or deactivated for the corresponding channel by pressing any numeric key while the cursor is located on the ' $\lambda 3$ ' or ' $\lambda 4$ ' field.
- 4) **Ratio plot 2** (the symbol λ_1^2 is shown above the decimal point of the absorption signal) The display shows the result of the division of signal (channel 2) by signal (channel 1).
Ratio plots can be selected on channel 3 and 4. They are activated, changed or deactivated for the corresponding channel by pressing any numeric key while the cursor is located on the ' $\lambda 3$ ' or ' $\lambda 4$ ' field.

Fig. 22 Signal negation and ratio plots: Setting output options and indication on the display.

Ratio plots contain information on **peak purity**. It is calculated by dividing the intensity values at two different wavelengths. It is possible to measure the absorption at one wavelength and to simultaneously monitor the ratio of absorption between this wavelength and another one. Ratio plots can be accessed at either analog out connector or via digital data acquisition.

The scaling of the analog outputs can be adjusted to different values of **Absorption Units** or to different voltage values. For details see section 'Analog out' on page 18 in the 'SETUP Menu' chapter.

EVENTS

The EVENT outputs (Event 1, Event 2 and Valve) which belong to the remote control sockets on the rear side of the instrument provide time programmed signals to other instruments. All features of the remote control sockets are described in the chapter 'Connecting other Instruments to the Smartline Detector 2600' on page 39.

When working without a time program, the events can be controlled manually or be used for fraction collection (see page 32).

Setting event outputs is done by entering numbers 0, 1 or 2 (0 = not active, 1 = active, 2 = pulse) in the event fields of the main menu. Here position one stands for Event 1, position two for Event 2 and position three for Valve.

Programming the Smartline Detector 2600

Creating Programs

Time displays are given in minutes with decimal calculation of seconds, e.g. 0.30 min = 18 s. Time resolution is 1/100 min. All programs are saved permanently and are present even after switching the instrument off.

Wavelengths and events can be changed as a function of time. Setup parameters such as scaling of the analog outputs cannot be time programmed. All options for data output like ratio plots (see page 28) can be defined for entire programs and will be used when loading a program.

The Smartline Detector 2600 can store up to twenty programs. By moving the cursor to the VIEW menu, you can check the availability and allocation of program numbers.

Use the EDIT mode for modifying or creating programs. Move the cursor to the **Edit Prog** field and enter the desired program number with numeric keys. Any numeric values can be overwritten. Pressing any arrow key confirms the entered value.

The first line of a program is the start time at 0 min. This time cannot be changed. Message **Start time is fixed** appears on editing this line. The maximum value for time is 999.99 min.

To edit a new time line, move the cursor to the **Time** field and press the ∇ key. If this line is empty, the **Time** field will show: '***.**'.


Fig. 23 Empty line in program no. 3

You can delete a program line if you enter '000.00' in the **Time** field and press any arrow key. The message **Delete this line? Confirm by cursor** will be displayed. You can confirm deletion if you press any arrow key. If you press a numeric key, or wait for 60 seconds, deletion will be aborted.

SOP 4 Time Programming of Wavelengths

Use this instruction for programming the absorption wavelength(s) as a function of time:

- 1) The cursor is in the main menu or it is moved back to the main menu.
- 2) Move the cursor to the **Edit Prog** field and enter a number for the new program (here: 03). The cursor goes automatically to time 000.00 which is the fixed start time (see **section A** in Fig. 24).
- 3) Press \triangleright to go to field $\lambda 1$. Enter the desired wavelength (e.g. 220 nm). Move the cursor to field $\lambda 2$ and enter the second desired wavelength (e.g. 254 nm). The same is done for channel three and four ($\lambda 3$ and $\lambda 4$). Entering nothing or 000 in a wavelength field deactivates this channel and the **Signal** field shows '---.'. The display should now resemble **section B**.


Fig. 24 Five steps to create a program

- 4) Press ∇ to create a new program line. The display should resemble **section C** with a blinking cursor on the first '*' of the **Time** field.
- 5) Now press 1 and 0 and then \triangleright to define time step 10.00 min. Enter 240 nm in field $\lambda 1$ and confirm with ∇ . The value for $\lambda 2$ (254 nm) and all others are kept constant as long as nothing else is entered. The display should now resemble **section D**.
- 6) Repeat the last two steps to define time step 20.0 min and enter 280 nm for $\lambda 1$. The display should resemble **section E**.

To run a program it must first be loaded. (See chapter 'Program Execution' on page 35.)

Event Programming

In addition to four detection channels two event outputs (one relay and one open collector) and one valve output (for details see page 39) can be programmed. Available settings are '0' (output deactivated), '1' (output activated) and '2' (pulse, duration: 500 ms).

During program execution, programmed events are visible only when they are active.


Fig. 25 Event programming

Fraction Collector Control

By using the fraction events (FRACTION EVENT1,2) and the fraction valve output (FRACTION VALVE) a simple, software independent fraction collector control can be implemented. One event controls the changing of fractions (this event is active only when a certain signal level is exceeded – peak detection!) and in addition a selection valve can be switched between *sample* and *waste* by the fraction valve output (this output switches the valve to *sample* when a certain signal level is exceeded and back to *waste* when the signal drops below this level).


Fig. 26 Submenus for fraction collection control

In the **DELAY** field a delay time (in seconds) may be entered that corresponds to the period of time a substance needs to get from the detector's flow cell to the fraction collector or fraction valve. This period depends on the volume between cell and fraction collector and on the flow rate.

Fraction events and fraction valve output are active only when they are switched **ON** from their respective submenus (Fig. 26) and this event must also be set to '1' or '2' (pulse) in the respective **Event** field of the main menu (logic AND operation). If a fraction event is set to '**OFF**' in the submenu or if a '0' is entered in an **Event** field, the fraction collection control is deactivated. *Time slots* can be created in programs by setting an event to '1' or '2' in the main menu for only certain periods of time.

Instead of a fraction collector, the KNAUER 7-Port 1-Channel (1/16" and 1/8"), 13-Port 1-Channel (1/8") or 17-Port 1-Channel (1/16") switching valves can be used for very flexible fraction collection. The control of these valves is identical to the control of a fraction collector.


Fig. 27 Necessary connections using KNAUER switching valves for fraction collection

If two switching valves are combined one valve (7-, 13- or 17-port 1-channel) can be used for the changing of fractions and the other one (6-port 3-channel or 7-port 1-channel) switches between *waste* and *sample* (Recycling Valve). The necessary electrical connections are shown in Fig. 27. As there are three event outputs available even more complex combinations of switching valves are possible.

SOP 5 Event Programming for Fraction Collector Control

Use this instruction to control a fraction collector via the Smartline Detector 2600 EVENT socket using signal level **and** a defined time window:

- 1) Enter the SETUP menu and proceed to the submenu **FRACTION LEVEL** (see page 19).
- 2) Enter the desired value for the signal level, e.g. 0,032 AU, from which on the event outputs should be activated.
- 3) Enable fraction event 1 by changing from 'OFF' to 'ON' in the **FRACTION EVENT1** submenu. Afterwards the desired DELAY can be entered if necessary.
- 4) In case an additional fraction valve is used the fraction valve output is enabled by changing from 'OFF' to 'ON' in the **FRACTION VALVE** submenu. Afterwards the desired DELAY can be entered if necessary.

- 5) Go back to the main menu.
- 6) Move the cursor to the **Edit Prog** field and enter a free program number (here: 07). Press ∇ to confirm.
- 7) Now program number 07 can be created. For this the desired wavelengths (here: 254 nm and 280 nm) for time 000.00 min are entered. The display should now resemble **section A** in Fig. 28.
- 8) After pressing ∇ the next program line appears. There the time '005.00 min' and the events '201' (activates fraction event 1 and fraction valve, which means that from 5 min on these two outputs are active as soon as the minimum signal intensity (fraction level) is exceeded) are entered. The display should resemble **section B**.
- 9) Press ∇ until the next program line is reached. Enter 010.00 in the **Time** field and '000' in the **Event** field. Consequently the events are deactivated after 10 min. The display should now resemble **section C**.

Event 1 (relay output) and the fraction valve output are activated between 5.00 min and 10.00 min. Every peak that exceeds the signal value which is set in the SETUP menu (submenu FRACTION LEVEL) (in this case 0.032 AU) within this time slot will produce a signal at event output 1 (pulse that advances the fraction collector to the next position) and at the fraction valve output (signal goes high to switch a valve from *waste* to *sample*).

Peaks that exceed the minimum signal value before 5 min or after 10 min do not produce any signal at any output.


Fig. 28 Program lines for program no. 7 (example)

The value for FRACTION LEVEL chosen in the SETUP menu, will always apply for all programs, it cannot be set individually for single programs.

Program execution

Programs must be activated prior to execution. Enter the desired program number in the field **Load Prog** and confirm by pressing an arrow key. The instrument's settings are adjusted to the programmed settings. An autozero is performed automatically every time a new program is loaded.

A program is started by moving the cursor to the **ON** field in the lower left corner of the display. Pressing '1' then starts the program. The display switches to **RUN**, and the elapsed time is displayed in the **Time** field. Program execution can be interrupted by pressing '0'. It is stopped with the actual values and the message **HOLD** is displayed.

Program execution can either be continued (return to **RUN**) by pressing '1' or aborted by pressing '0'. When a program is aborted the instrument returns to the start settings (first program line). These settings are also activated after a complete run (program) is finished.

A program can also be started via the **START IN** input at the instrument's rear panel. In case a trigger signal is applied to this input, the system switches into **RUN** mode.

Creating Program Links

Program links are cross references to existing programs, which can be defined and edited similarly to programs themselves. Links cannot reference other links.

When editing, mode links are defined and loaded as the reserved program numbers 21 to 29, where program 21= L1, 22 = L2 etc.

Links have line numbers from **No01** up to **No99**, every line containing one program, the number of repetitions and a programmable wait-status. The Smartline Detector 2600 can store up to 9 different links, but the maximum number of programs used in all links is limited to 100. Each program may appear in any frequency in any link.

Programming Links

Enter the number of the desired link in the field **Edit Prog**, for example 22 for link L2 and press any arrow key to confirm. **Part A** in Fig. 29 is showing that step, the cursor blinking at **No****. Now press any numeric key to generate a new, continuously numbered line. In this case **No01** will appear.

Move the cursor to the **P** field by pressing \blacktriangleright and choose the desired program number, e.g. 03. Move on to the **R** field to put in the number of repetitions for the selected program (e.g. 05). Finally, in the **Wait** field the values **1** or **0** can be entered. For Wait = 1 the detector will wait for an external start signal or for an operator pressing **1** with the cursor on the **RUN** field until this line will be performed. During waiting the instrument is switched to the wait-status. If Wait is set to **0**, the program lines of the link are running continuously without any break. This step is shown in part B of Fig. 29. In this example program 03 would run five times and then the next line of the link program would follow.

After creating a new line number a new link program line is made. Here the execution of another program may be added to the first program.

Between the lines of the LINK program you can move with the \blacktriangle and \blacktriangledown keys. From the first line and from the last line the main menu can be reached by pressing the \blacktriangle key or the \blacktriangledown key respectively.


Fig. 29 Editing program links

Running Program Links

Before a program link can be executed, it has to be loaded. For this purpose enter the desired link number in the **Load Prog** field, e.g. 22 for L2 and confirm by pressing any arrow key. A program link is started exactly like a program. A program is started by moving the cursor to the **ON** field in the lower left corner of the display. Pressing '1' then starts the program. The display switches to **RUN**, and the elapsed time is displayed in the **Time** field.

During the execution of a link the LINK menu is available in addition to the common display. In the LINK menu the momentary status of the running link is shown. Move cursor to the **Load Prog** field and press \triangleleft to enter this menu page. For more explanations refer to 'LINK menu' on page 23. To leave the LINK menu press \triangleright .

Deleting Programs and Program Links

To delete a program or a link, it has to be loaded first (enter the program number in the **Edit Prog** field and confirm with any arrow key). Entering 0 in the **Edit Prog** field and pressing an arrow key will produce the message **Delete this program? Confirm by cursor** in the display. Pressing any arrow key will delete the loaded program or link. By pressing any numeric key or waiting 60 seconds this operation will be cancelled (the program or link will not be deleted).

Wake Up Program

Program no. 30 is reserved for the Wake Up function, enabling the start of the lamp and/or a wavelength program at a preset time or date. The Wake Up function can be used e.g. for ensuring switch on and preheating of the lamp after switching it off overnight. An example for the Wake Up program is shown below.


Fig. 30 Programming the Wake Up Program

The number of the program and the time and date for wake up can be entered with the numeric keys. In the **DELAY** field a period of time (up to 99 minutes) may be entered that the instrument will wait after the lamp has been started and before the program is started.

Using the Wake Up Program

On loading the Wake Up program the question: 'go to sleep mode?' will be displayed (**section A**). Pressing any arrow key will then switch the lamp off and start the Wake Up program. The display shows: 'SLEEP MODE', as well as the current time and the time for wake up (**section B**).


Fig. 31 Display showing the Wake Up program

As soon as the preset time is reached, the lamp is switched on, the wavelength set according to the chosen program and an autozero is performed (the chosen program is started). If a delay was set, the program will be started later.

Scanning UV Spectra

The Scan Menu

The scan menu is accessed upon pressing the **SCAN** key. This menu enables choosing the scan number (1-4) to save the scan data and serves for giving out the stored scan data. The wavelength range is preset in the hardware and can not be changed.

A	SCAN no.:1	OUT	B	SCAN no.:1	OUT
	◆ range	192. . . 510		◆ range	430. . . 710

Fig. 32 The Scan Menu (A deuterium lamp, B tungsten halogen lamp)

The Smartline Detector 2600 is capable of storing up to four different scans in the internal memory. The four scans are numbered from **no.:1** to **no.:4**. Once saved in the memory, they can be activated and given out via the analog outputs e.g. to a chart recorder.

The number of the scan to be performed or given out is selected by moving the cursor with **▷** to the **SCAN** number field, e.g. **no.:1**. The number can be changed in this position with **▽** or **△**.

Performing Wavelength Scans

For acquiring scan data, enter the scan menu by pressing the **SCAN** key once. When pressing the **SCAN** key a second time, a scan is performed immediately and stored under the preset scan number. As a confirmation the display shows e.g. '**Scan No. 1 acquired**'.

Spectra can be recorded in standalone mode. That means also without KNAUER chromatography software. Spectra scans can not be programmed (in standalone mode).

For taking different scans during a run, please be sure to adjust the scan number between scans. Stored scans are overwritten without any message, when the **SCAN** key is pressed twice.

Stored spectra can only be erased by overwriting them with new spectra.

Output of Scan Data

The stored scan data can be given out via **analog output 1** and with defined data rates. Time scaling of the data can be adjusted between 1 and 50 steps per second (that is approximately 1-50 nm/s) in the **SETUP** menu of the detector, see section **SCAN OUT** on page 18.

SCAN OUT
◆ speed: 1 nm/s

Low output speed settings lead to enhanced spectral resolution of the data output.

To release the output of scanned data, move the cursor using the **▷** key to the **OUT** field in the scan menu and press the **△** or the **▽** key to start the scan output. Make sure that the desired scan number is chosen. Since every scan is kept in internal storage it can be edited or given out several times. Scans are deleted by overwriting with new scan data. During scan data output the recent wavelength is shown on the display:

SCAN OUTPUT 198 nm

Connecting other Instruments to the Detector 2600

Using the Remote Sockets

The remote sockets are located on the rear panel of the Smartline Detector 2600, pos. **3** and **5** in Fig. 7. They serve to send or/and receive signals to/from other instruments. For example the start signal from an injection valve or an autosampler can be put on the detector's **START** input. The output **Event 1** (relay) can for example be used in combination with the detector's fraction collection control functionality. A KNAUER 16-way switching valve or a fraction collector can be connected to this output and is advanced by one position each time a certain signal level is exceeded (see page 32).

Please avoid touching the electrical contacts of the socket lines. Electrostatic discharges when touching the contacts could damage the electronics of the device.

Connecting the Remote Sockets

The remote sockets offer two Event outputs one Error input, Error output, Autozero input, Start input, Start output and one Fraction valve output.

Remote Socket 1


Fig. 33 Remote socket 1

- EVENT 1 Relay:**
Possible settings:
contact open (0)
contact closed (1)
pulse (contact is closed for 500 ms) (2)
EVENT 1 output is switched from programs and/or from the main menu.
- EVENT 2 Open Collector (max. 100 mA) :**
Possible settings:
low (no current) (0)
high (current up to 100 mA) (1)
pulse, (current up to 100 mA for 500 ms) (2)
EVENT 2 output is switched from programs and/or from the main menu.
- GND Ground for EVENT 2**
- ERROR IN** A short circuit (via relay or Open Collector) to **GND** causes an error message on a connected PC.
- ERROR OUT Open Collector.** This output is active as long as an error is shown on the display – e.g. when the lamp does not start (Error: 0 V, no error: +5 V).
- ERROR GND Ground** for ERROR IN and ERROR OUT

Remote Socket 2


Fig. 34 Remote socket 2

AUTOZERO A short circuit to **GND** causes an autozero for all channels.

GND **Ground** for **AUTOZERO**

START IN A short circuit to **GND** causes the start of a program (in case a program is loaded) . **START IN** can also send a signal to the KNAUER Software .

START OUT **Open Collector**. This output stays active (+5 V) for 500 ms when a program or link is started.

START GND **Ground** for **START IN** and **START OUT**

VALVE + **+24 V** for **VALVE -**

CAUTION: This output is used only for connecting 24 V valves that have no own power supply. KNAUER switching valves are connected to GND and Valve – (see Fig. 27).

VALVE - **Open Collector (max. 200 mA / 24 V) :**

Possible settings:

low (no current) (0)

high (current up to 200 mA) (1)

pulse, (current up to 200 mA for 500 ms) (2)

The **VALVE** output is switched by a program or via the main menu.

Assembling WAGO Plugs

For the electrical connections to other instruments via the remote control strip enclosed 8-fold WAGO plug strips are used. They may be mounted as follows.


Fig. 35 Assembling plus

SOP 6 Assembling WAGO plugs

1. Insert the rounded end of the lever latch into the square opening of the selected connector of the plug.
2. Press the lever down as indicated by the arrow in Fig. 35.
3. Insert the non insulated end of the cable into the opening under the lever.
4. Release the lever and remove the lever latch from the plug.

The cable is now firmly anchored in the plug.

The Integrator Outputs

Fig. 36 Integrator output (1-4)

The Smartline Detector 2600 has four analog integrator outputs for up to four detection channels. At an integrator output the recent signal value is represented by an analog voltage. In the SETUP menu, submenu ANALOG OUT FULL SCALE, the signal can be scaled in 16 steps (see page 18). Default setting is 1 AU = 1 V. The integrator outputs' maximum output voltage is $\pm 2,5$ V (equals $\pm 2,5$ AU at default setting). If necessary, an individual offset voltage (max. 999.99 mV) can be selected for every channel from the SETUP menu, submenu ANALOG OUT OFFSETS. The outputs are connected via the enclosed analog cables (cinch-cinch) or via optional cables with other connectors.

Output 1 is also used to send spectra that were stored in the instrument (see spectra output on page 38).

All integrator outputs are active in standalone mode as well as under control of KNAUER Software.

RS-232 Serial Interface

Fig. 37 RS-232 Interface

The RS 232 serial interface on the rear side of the device enables digital data transfer between the Smartline Detector 2600 and a PC equipped with e.g. KNAUER ChromGate® (from Version 3.1), ClarityChrom®, or EuroChrom® (from Version 3.05) software. Detector and PC are connected via the enclosed RS 232 female/female cable.

Alternatively, the instrument may be controlled using a simple ASCII protocol.

Contrary to other KNAUER instruments, this detector's RS 232 interface runs at 115 kBaud (rather than at 9.6 kBaud). This is necessary due to the high amount of data that is being sent.

The KNAUER Net protocol is not supported.

USB Interface


Fig. 38 USB Interface

The USB interface is not functional and may be implemented at a future date.

Software Control of the Detector 2600

The full capabilities of the Smartline Detector 2600 within an HPLC system are accessed only under operation with one of the KNAUER software packages ChromGate® (from Version 3.1), ClarityChrom®, or EuroChrom® (from Version 3.05)®. When controlled by the software the keypad is locked and the display just shows: **'REMOTE CONTROL'**.


Fig. 39 ChromGate® HPLC software

In combination with the software, all features of the diode array technology on which this detector is based can be used. Depending on the software, functions such as continuous spectra scanning, 2D chromatograms, 3D chromatograms, peak purity functionality, spectral library functionality and more become available.

For KNAUER ChromGate® Software, an additional *scan license* is required.

For more detailed information about the software, its features and how to work with it, please have a look into the particular software manual.

Simple Maintenance

Especially for routine checks of the instruments you may use the KNAUER OQ-documents (available from your local dealer). An automated OQ-check is available with KNAUER ChromGate® Software (from V3.1 on).

Control of the Lamp functionality

The deuterium lamp used with the Smartline Detector 2600 has an extended life time to ensure long-time functionality and reliable measurements with low noise and baseline drift as well as high sensitivity. The actual working time of the lamp depends on different factors, like the number of lamp starts, the average burning time and the requirements concerning noise and sensitivity.

To check the functionality of the lamp, the intensity value **sig1** (found in the SETUP menu, submenu SIGNAL) gives useful information but only if the wavelength is set to 240 nm for the deuterium lamp or 520 nm for the tungsten halogen lamp and if the instrument is operated with a dummy cell. The **sig1** value equals the light intensity which gets to the diode array and can be used to estimate the lamp quality.

SIGNAL:	sig1:0.4544
◆ absorption	sig2:0.0000

Fig. 40 SIGNAL menu

On delivery of the detector, **sig1** is set to the range of 0.7 to 0.9 measured with the dummy cell at a wavelength of 240 nm or 520 nm respectively.

We recommend to check the sig1 value at regular intervals under the conditions mentioned above (dummy cell, $\lambda = 240$ nm or 520 nm respectively, integration time < 40 ms). This applies especially in case higher noise levels or decreased sensitivity are observed on working with the photometer K-2600. If these observations coincide with a **sig1** value of approx. 0.1 or less, a new deuterium lamp should be installed (see 'Changing the lamp' on page 45).

Checking Wavelength Accuracy

The KNAUER Smartline Detector 2600 is equipped with an integrated holmium oxide filter. This enables manual or automated wavelength accuracy check.

Holmium oxide glass has an absorption spectrum with characteristic absorption maxima. By comparing the positions of the detected maxima with the values in the enclosed holmium filter certificate the detector's wavelength accuracy can be determined.

For checking wavelength accuracy of the Detector 2600 according to SOP 7 a spectrum is recorded with the integrated holmium oxide filter moved into the optical path.

Alternatively the automatic wavelength accuracy check of the OQ procedure in KNAUER ChromGate® Software can be used.

Due to the completely different algorithms and recording conditions the results as well as the limits of the wavelength accuracy checks (standalone or software controlled) differ significantly.

SOP 7 Wavelength accuracy check

1. Install the supplied dummy cell (test cell).
2. Switch on the detector and wait until the deuterium lamp has started.
3. Move to the GLP menu, submenu **Wavelength validation** and press the \triangleright key. The following menu appears:

Perform validation?
cursor – YES, other - NO

4. After confirming with one of the arrow keys (cursor) the automatic validation will be run. The display shows:

Performing validation
please wait...

5. After successfully checking four defined holmium lines (that means all deviations were $\leq 1\text{nm}$) the result is displayed and after pressing any key the instrument is back in the submenu Wavelength validation.

Validation OK
Press any key

6. By pressing the ∇ key the results of the last wavelength validation are displayed in detail.

HO λ	Meas. λ	
◆ 241.5	240.8	1 OK

HO λ	Meas. λ	
◆ 240.9	240.6	1 OK

HO λ	Meas. λ	
◆ 279.2	278.8	2 OK

HO λ	Meas. λ	
◆ 445.7	445.3	2 OK

HO λ	Meas. λ	
◆ 361.0	360.6	3 OK

HO λ	Meas. λ	
◆ 536.5	536.1	3 OK

HO λ	Meas. λ	
◆ 446.0	445.4	4 OK

HO λ	Meas. λ	
◆ 637.7	637.2	4 OK

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tungsten halogen lamp

For four holmium lines the defined wavelength (**HO λ**), the measured wavelength (**Meas. λ**) and the result (**OK** or **Failed**) are shown.

In case the validation fails repeatedly, the instrument should be recalibrated.

Changing the lamp

The exchange of a deuterium lamp by a tungsten halogen lamp or vice versa is not possible.

Remove the power plug before opening. Please let the lamp cool down for at least 15 minutes after switching it off.

Do not touch the glass of the lamp. Should you touch it accidentally, clean it thoroughly with a lint free cloth and isopropanol.

Fig. 41 Deuterium lamp with plug

SOP 8 Changing the deuterium lamp

1. Open the detector's front door.
2. Unscrew screw (1) with the supplied Allen key (size 3).


Fig. 42 Detector with open door

3. Pull the optical bench out of the housing.

Fig. 43 Detector with optical bench

4. Pull the 3-pole plug (2) of the old lamp. Unscrew the two screws in the lamp socket (3) with a screwdriver and remove the whole lamp including its cables from the instrument.


Fig. 44 *Optical bench with lamp*


Fig. 45 *Optical bench, lamp removed*

5. Insert the new lamp assuring that it is correctly seated in the guiding slot.
6. Screw it in securely and connect the 3-pole plug into the 3-pole socket (2).
7. Increase the lamp counter in the GLP menu, submenu LAMP by one. The operating time counter will be set to zero automatically.
8. Put the optical bench back into the instrument carefully.
9. Tighten the screw (1).
10. Switch the detector on again and check the intensity values in the SETUP menu, submenu SIGNAL.

When a new lamp is installed, it takes about 24 hours to reach the optimum working conditions.

SOP 9 Changing the tungsten halogen lamp

Fig. 46 Tungsten halogen lamp

1. Open the detector's front door.
2. Unscrew screw (1) with the supplied Allen key (size 3).
3. Pull the optical bench out of the housing.
4. Pull the 2-pole plug of the old lamp. Unscrew the two screws in the lamp socket with a screw driver and remove the whole lamp from the instrument including its socket and cables.
5. Loosen the pin screw (4) and pull the lamp out of its socket.
6. Insert the new lamp assuring that it is correctly seated in the guiding slot.
7. Screw it in securely and connect the 3-pole plug into the 3-pole socket (2).
8. Increase the lamp counter in the GLP menu, submenu LAMP by one. The operating time counter will be set to zero automatically.
9. Put the optical bench back into the instrument carefully.
10. Tighten the screw (1).
11. Switch the detector on again and check the intensity values in the SETUP menu, submenu SIGNAL.

Cleaning the flow cell

Noisy baselines and low sensitivities may be due to a dirty flow cell. This may also be indicated by a low value for **sig1** (240 nm or 520 nm respectively) in the SETUP menu, submenu SIGNAL when flushing the cell with pure solvent. In most cases it is sufficient to purge the flow cell according to the following SOP.

SOP 10 Purging the Flow Cell

1. Purge the flow cell using one of following solvents: sodium dodecyl sulfate (SDS), 1M HCl, 1M NaOH, ethanol, or acetone.
2. Run the solvent through the flow cell using a syringe and leave it in the cell for approximately 5 minutes.
3. Rinse extensively with water and then blow dry using a gentle stream of pure nitrogen.

Never dry with compressed air from a "house" line as this will contain microdroplets of oil that will coat the cell.

When the cell is not in use, disconnect the flow cell and clean out traces of salt and protein with a syringe filled with distilled water.

Before storing the flow cell, inject a dilute solution (10-25%) of ethanol or isopropanol to prevent microbial growth.

In case flow cell purging does not provide sufficient success, all flow cells can easily be disassembled for cleaning the lenses for example in an ultrasonic bath.

Analytical flow cells


Fig. 47 Sectional view of an analytical flow cell

SOP 11 Disassembling and Cleaning of an Analytical 10 mm Flow Cell

1. Unscrew the outer threads with the 3 mm hexagonal spanner enclosed in the flow cell's delivery.
2. Remove the black gasket that carries the lenses with a pair of tweezers or by gently tapping it on a clean surface. The lens is embedded in the gasket and sealed against the flow path with a PFA seal. This seal should be changed every time when disassembling the flow cell.
3. Take out the lenses and clean them by wiping them with a soft cloth or with an appropriate solvent in an ultrasonic bath. Be careful not to touch the clean lenses directly.
4. Reassemble the cell in the reverse manner, making sure that the PFA seal does not block the light path.
5. Tighten the outer threads carefully with the spanner in order not to damage the lenses.

Analytical 3 mm Flow Cells


Fig. 48 Sectional view of an analytical 3 mm flow cell

SOP 12 Disassembling and Cleaning of an Analytical 3 mm Flow Cell

The 3 mm flow cells (A4042, A4045) have a rod shaped light guide and one concave lens.

1. Unscrew the outer thread with a hexagonal spanner.
2. Take out the stainless steel shutter.
3. Get hold of the gasket holder with the light guide by grasping it with a pair of tweezers, using the indentations on the outer side of the gasket holder.
4. Push out the light guide and strip the PTFE sealing ring in order to clean the lens.
5. Unscrew the second outer thread with the 3 mm hexagonal spanner enclosed in the flow cell's delivery.
6. Remove the black gasket that carries the lenses with a pair of tweezers or by gently tapping it on a clean surface.
7. Remove the PEEK plate.
8. Get hold of the gasket holder with the light guide by grasping it with a pair of tweezers, using the indentations on the outer side of the gasket holder.
9. Push out the light guide and strip the PTFE sealing ring in order to clean the lens.
10. Reassemble the cell in reverse order. Use a new PTFE sealing ring after every disassembly to ensure the consistence of the flow cell.

Preparative Flow Cells


Fig. 49 Sectional view of a preparative flow cell

SOP 13 Disassembling and Cleaning of a Preparative Flow Cell

The preparative flow cells (A4066, A4067, A4069, A4095) have a rod shaped light guide instead of the concave lens in the analytical cells.

1. Unscrew the outer thread with a hexagonal spanner.
2. Take out the stainless steel shutter and the PEEK distance plate (not present in A4069, A4095).

3. Get hold of the gasket holder with the light guide by grasping it with a pair of tweezers, using the indentations on the outer side of the gasket holder.
4. Push out the light guide and strip the PTFE sealing ring in order to clean the lens.
5. Reassemble the cell in reverse order. Use a new PTFE sealing ring after every disassembly to ensure the consistence of the flow cell.

Adjusting the Path Length of Preparative Flow Cells

SOP 14 Changing the Path Length of Preparative Flow Cells

Path lengths of preparative flow cells (A4066, A4067 - 1/8" only) can be adjusted to 2, 1.25 and 0.5 mm. On delivery, the path length is set to 2 mm. To reduce the path length to 1.25 or 0.5 mm, follow the instructions given:

1. Unscrew the outer thread with the 3 mm hexagonal spanner.
2. Take out the stainless steel shutter and the PEEK distance plate.
3. Remove the distance plate, put back the stainless steel shutter and refasten the thread carefully.

With the distance plate removed, the rod shaped light guide is pushed further into the flow cell (0.75 mm), thus resulting in a reduced path length of 1.25 mm. To reduce the path length further to 0.5 mm, follow the same procedure on the other side of the cell.

To extend the path length again in steps of 0.75 mm, the PEEK distance plates have to be inserted again:

1. Loosen the outer thread, remove the stainless steel shutter and take out the gasket holder by using a pair of tweezers.
2. Push the light guide approximately 1 mm to the outside to enlarge the path length. Use a clean cloth and do not touch the light guide with the fingers.
3. Put the gasket holder back into the cell.
4. Insert the PEEK distance plate again and then the shutter.
5. Fasten the outer thread carefully.

When fastening the threads, the rod shaped light guide is pushed back in the correct position inside the cell. Inserting a distance plate thus enlarges the path length for 0.75 mm. It is not necessary to change the PTFE sealing ring when adjusting the path length.

Error Messages

Error Messages, Descriptions and Codes

Error message	Error description	Error code
Shutter wheel error	Shutter test not successful	1
Lamp intensity too low	Lamp did not start or light intensity too low	2
λ validation error	Wavelength validation with internal holmium oxide filter failed (deviation < ± 1 nm)	3
External error	A signal was detected at the Error Input	4
Communication error	Error in communication between detector and PC	5

Table 5 Error messages

Error Messages and their Effects

Error message	Reaction of the instrument	Error code saved (GLP)	Error Out active
Shutter wheel error	Fatal error. Instrument must be switched off.	Yes	Yes
Lamp intensity too low	Fatal error. Instrument must be switched off.	Yes	Yes
λ validation error	Press any key to continue. The error message will be repeated from time to time.	Yes	Yes
External error	Press any key to continue. The error message will be repeated from time to time.	No	Yes
Communication error	The instrument tries to reinstall communication. Otherwise the error message will be shown.	Yes	Yes

Table 6 Error messages and their effects

Error Messages and their Reasons

Error message	Possible reason	Correction
Shutter wheel error	Optical bench is not completely pushed into the instrument.	Please slide the optical bench completely into the instrument and screw it tight.
	Mechanical, electrical or electronic problem	Please call KNAUER service department.
Lamp intensity too low	Lamp lifetime exceeded	Please replace the deuterium lamp.
	Flow cell dirty	Please clean the flow cell.
	Problem in the optical system	Please call KNAUER service department.
λ validation error	Not enough light	Run the automatic integration time setting with the dummy cell installed.
		Replace the lamp if sig1 << 0.1 (at 240 nm or 520 nm, with dummy cell).
	Too much light	Run the automatic integration time setting with dummy cell.
	Instrument is decalibrated	Please call KNAUER service department.
External error	Error outside of the detector	Please remove the error.
Communication error	Communication cable loose or defective	Please check the communication cable. Replace the cable if necessary.
	Software error	Please restart the software and if necessary also the PC
	Error of the instrument	Please call KNAUER service department.

Table 7 Error messages and their reasons

Other Problems and possible Causes

Problem	Possible causes	Solution
Low sig1 value at 240 nm or 520 nm also with dummy cell (SETUP menu, submenu SIGNAL)	Integration time too low	Increase the integration time or run the automatic integration time setting.
	Old lamp	Please replace the deuterium lamp.
Low sig1 value at 240 nm or 520 nm with flow cell, <u>not</u> with dummy cell (SETUP menu, submenu SIGNAL)	Cell is dirty	Please clean the flow cell.
	Eluent absorption too high	Please replace the eluent. Please try to work at a different wavelength (less absorption by the eluent).
Spikes in the baseline	Bubbles	Please check the tubes on the intake side of the pump for leakage.
		Please use a degasser.
		Please try to purge the cell at high flow rate.
		Please check with dummy cell.
High drift	Cell is leaking	Please remove the leakage.
	Heavy temperature changes	Please work at constant room temperature..

Table 8 Other problems

Spare Parts and Accessories

Flow cells for the Smartline Detector 2600

Analytical Flow Cells

Order No. Cell type	Path Length (mm); Connector	ID channel (mm)	Volume (μ l)	Material	Flow Range (ml/min)	Maximum Pressure (bar)
A4061	10 mm; 1/16 "	1.1	10	stainless steel, with heat exchanger	20	300
A4042	3 mm; 1/16 "	1.0	2	stainless steel	50	300
A40 45	3 mm; 1/16 "	1.0	2	PEEK	50	30

Preparative Flow Cells

A4066	0.5/1.25/2 mm 1/8"		1.7/4.3/ 6.8 μ l	stainless steel	1.000	200
A4067	0.5/1.25/2 mm 1/8"		1.7/4.3/ 6.8 μ l	PEEK	1.000	100
A4068	0.5/1.25/2 mm 1/4"		1.7/4.3/ 6.8 μ l	stainless steel	10.000	200
A4069	0.5 mm 1/16"			stainless steel	250	200
A4095	0,5 mm 1/16"			PEEK	250	100

U-Z View™ Micro Flow Cells

A4091	8 mm 1/16"	0.150	0.140	fused silica	0.10	500
A4092	8 mm 280 μ m	0.015	0.035	fused silica	0.01	500

Spare Parts and Accessories, Order Numbers

A4071	Deuterium lamp
A4072	Tungsten halogen lamp
A4073	Tungsten halogen lamp including socket
M1479	Power supply cable
A0895	RS-232 connection cable (9 pin, female/female)
M0205	WAGO plug strip (8 x)
M0156	WAGO lever latch
M1588	Analog connector cable
G1023	Integrator cable
A1467	10 pin ribbon cable
A1131	Repair kit for analytical flow cells
A1132	Repair kit for preparative flow cells
A4124	Test cell (dummy)

Specifications

Lamp	Deuterium or Halogen
Wavelength range ^{*)}	190 - 510 nm (D ₂) 430 - 710 nm (tungsten halogen)
Detection type	Diode array
Dot pitch	1.25 nm
Diode array	256 Diodes
Bandwidth	4-25 nm (selectable)
Acquisition channels (with software / standalone)	4 / 4
Analog outputs (scalable)	4 x ± 2,5 V, 20 bit, scalable in 16 steps
Spectra to be stored	4
Wavelength accuracy	≤ 1 nm
Wavelength precision	≤ 0,5 nm
Wavelength validation	Automatic, with internal holmium oxide filter
Noise according to ASTM E1657-94, optical path 10 mm, 1 ml/min methanol, λ = 254 nm, Δλ = 4 nm (± 2 nm)	≤ 1·10 ⁻⁵ AU (D ₂), 2·10 ⁻⁵ AU (tungsten halogen)
Drift according to ASTM E1657-94	≤ 5·10 ⁻⁴ AU/h (D ₂), 3·10 ⁻⁵ AU/h (tungsten halogen)
Linearity	0 – 1.5 AU
Measurement range	0 – 2.2 AU
Time constants	0.1/ 0.2 / 0.5 / 1.0 / 2.0 / 5.0 / 10.0 s
Integration time range	13 – 200 ms
Control	RS 232 interface, remote connector, keypad
Programs	20 Programs, 9 Links, 1 Wake up
Display	LCD, 2 lines, 24 digits
Power Supply	90-260 V, 47 – 63 Hz, 70 VA
Dimensions (W x H x D)	226 x 135 x 390 mm
Weight	5.7 kg
GLP Support	Detailed trace report with operating hours of total, lamp; number of lamp ignitions, service information

^{*)} The range may be marginally shifted for individual instruments.

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Warranty statement

The warranty period of the Smartline Detector 2600 (UV or VIS) is 12 months beginning from the date of dispatch from Berlin. Operation inconsistent with manufacturer's instructions or damage caused by unauthorized service personnel are excluded from guarantee. Damage caused by blockages and wear and tear parts such as fuses and seals are not covered by the guarantee. Defective detectors should be sent to the manufacturer for repair.

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If we find a defect covered by the guarantee, repair or replacement, at our discretion, will be carried out free of charge. Packing and transport costs are borne by the purchaser.

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Declaration of conformity

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Manufacturer's name and address:
Herstellername und -adresse

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Smartline Detector 2600,

Order Numbers, Bestellnummern: **A 5200, A5201, A5210 and A 5211**

complies with the following requirements and product specifications:

- Low Voltage Ordinance (73/23/EWG)
EN 61010-1 (08/2002)
- Engineering Guidelines (89/392/EWG)
- EMC Ordinance (89/336/EWG)
EN 6100-3-2 (10/2006)
EN 61326-1 (10/2006)

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Index

- A**
- Accessories 54
 - Arrow keys 14
 - Auto zero 14
 - Autotest 13
- B**
- Bandwidth 18
- C**
- Capillary connection 26
 - ChromGate 42
 - OQ procedure 43
 - ClarityChrom 42
 - Cursor appearance 16
- D**
- Delivery 10
 - Deuterium lamp 45
 - Direct control 28
 - DYNASEAL 27
- E**
- Error Messages 51
 - Events 29
- F**
- Firmware 22
 - Flow cell 23, 54
 - Analytical 48
 - Cleaning 47
 - Preparative 49
 - Foil keys 14
 - Fraction
 - Event 19
 - Level 19
 - Valve 19
 - Fraction collector control 32
 - Front panel 11
- G**
- General description 8
 - GLP menu 20
- I**
- Integration time 19, 27
- L**
- Lamp 18
 - Changing 45
 - Functionality 43
 - LINK menu 23
- M**
- Main Menu 15
 - Maintenance 43
 - Menu structure 15
 - Micro flow cells 54
- O**
- Optical path 9
- P**
- Path length 50
 - Power supply 13, 28
 - Program
 - Creating 30
 - Deleting 36
 - Event 32
 - Execution 35
 - Links 35
- R**
- Rear Panel 12
 - Remote control sockets 39
 - RS-232 Interface 41
- S**
- Scan 14, 18
 - Scan Menu 38
 - SETUP Menu 17
 - Shutter 19
 - Signal 18
 - Software control 42
 - SOP 7
 - Spare parts 54
 - Specifications 55
- T**
- Time constant 17
 - Tungsten halogen lamp 47
- U**
- USB Interface 42
- V**
- VIEW menu 22
- W**
- WAGO-plugs 40
 - Wake Up Program 37
 - Warming up 28
 - Wavelength
 - Accuracy 43
 - Selection 28
 - Validation 21
 - Working time 21

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