K-7000

Vapor Pressure Osmometer

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User Manual / Handbuch

V7109, 03/2007

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Table of contents

| | |
|--|----|
| Design & Measurement Principle | 5 |
| Preparing the Vapor Pressure Osmometer K-7000 for Measurement..... | 6 |
| Front View of the Instrument | 6 |
| Electrical Connections & Rear Panel..... | 6 |
| Rear Panel Elements..... | 6 |
| Top View of the Instrument..... | 7 |
| Removing the Cell from the Instrument..... | 8 |
| Cleaning the Thermistor Probes..... | 8 |
| Insertion of Vapor Wicks..... | 8 |
| Insertion of the Measurement Cell..... | 8 |
| Operating your Vapor Pressure Osmometer K-7000..... | 9 |
| Switching the Vapor Pressure Osmometer K-7000 on..... | 9 |
| Using the Keys of the Vapor Pressure Osmometer K-7000 | 9 |
| Displays of the Vapor Pressure Osmometer K-7000..... | 10 |
| Main Menu | 10 |
| Selecting an appropriate Measurement Temperature..... | 10 |
| Using the Time Display..... | 11 |
| Identifying the Serial Number | 11 |
| Using the Calibration Mode | 11 |
| Using the Mode Actual Temperature..... | 12 |
| Resetting the Bridge Voltage..... | 12 |
| Practical Choice of Solvents..... | 12 |
| Specifications of the required Syringes | 14 |
| Performing Measurements | 14 |
| Using the syringes | 14 |
| Finding the Zero Point | 14 |
| Dropping Techniques..... | 14 |
| Measuring new solvent and new sample droplets..... | 15 |
| Exchanging the sample droplet only..... | 15 |
| Thermistor Adjustment | 16 |
| Required Stabilization Time..... | 16 |
| Three Different Measurement Modes | 16 |
| Determination of the Total Osmolality of Physiological Fluids..... | 16 |
| Performing a Calibration with NaCl of 400 mOsmol/kg..... | 18 |
| Evaluation of calibration data | 18 |
| Evaluating Measurement Data | 19 |
| Molecular Mass Determination for $M < 500$ g/mol..... | 19 |
| Molecular Mass Determination for $M > 500$ g/mol..... | 21 |
| Calibration..... | 21 |
| Measurement:..... | 21 |
| Operating the K-7000 with EuroOsmo [®] 7000 Software..... | 23 |
| Hardware Requirements and Installation of EuroOsmo [®] 7000 Software...23 | |
| Installation of the EuroOsmo [®] 7000 software from a CD ROM..... | 23 |
| Installation of the EuroOsmo [®] 7000 software from a floppy disk | 23 |
| Hardware installation | 23 |
| Explanation of the main screen | 24 |
| Software set up for the measurement | 27 |
| Performing measurements | 27 |
| Determination of the total osmolality | 28 |
| Determination of molecular mass < 500 | 28 |
| Determination of molecular mass > 500 | 29 |
| Data saving..... | 30 |
| Trouble Shooting Tips..... | 32 |
| Spare Parts & Accessories | 33 |
| Data Sheet of the Vapor Pressure Osmometer K-7000 | 34 |
| Warranty statement..... | 67 |
| Declaration of conformity | 68 |
| KNAUER Software License Agreement..... | 69 |

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Design & Measurement Principle

The Vapor Pressure Osmometer K-7000 has been designed to exactly measure the total osmolality of physiological fluids and to determine the average molecular mass of polymers in aqueous or organic solutions.

Two thermistors are part of a Wheatstone bridge. Differential measurements of currents due to temperature differences at the thermistors are accomplished. The thermistors are located in a cell where the gas phase is saturated with solvent vapor. The cell temperature is electronically controlled and maintained with an accuracy of $\pm 1 \times 10^{-3} \text{ }^\circ\text{C}$. One drop of the same pure solvent located on either thermistor represents a zero temperature difference or equilibrium of the measurement system.

The vapor pressure of any solution containing solutes is lower than the vapor pressure of a pure solvent. Hence, replacing one drop of pure solvent with one drop of a solution leads to a vapor pressure difference between the two droplets. This difference however is compensated as follows: Some vapor of the pure solvent that saturates the gas phase condenses on the droplet of the solution until the vapor pressures are balanced. The increasing vapor pressure of the solution droplet leads to an increase of temperature. Once equilibrium is reached a constant measurement value is achieved.

This particular ΔT between the thermistors is always proportional to the number of particles or number of moles dissolved in the solution. Thus either concentrations or the molecular mass can be determined, in the latter case however only if sample concentrations are known.

The measuring cell is placed in a thermostat consisting of an anodized aluminum block. The two thermistors immerse into a glass beaker enriched with vapor. The operating temperature can be selected in the range between 20 and 130°C. The cell is sealed with a circular cell cover and a teflon gasket.

The head of the osmometer has 6 directional channels where the syringes can be introduced to attach droplets to the thermistors. The syringes are introduced by slightly moving them downwards into individual channels. The needles of the syringes are equipped with springs. Direct contact with the probes should be avoided. The syringe ports are numbered from 1 to 4 for the solutions, the 2 ports without numbers are to provide both thermistors with pure solvent drops.

The sample ports are aligned with the left thermistor. The head of the osmometer will also be thermostatted. To preheat the syringes it is recommended to set the head temperature 2°C higher than in the measuring cell. This value of experience provides an excellent screening against outside temperature variation and a minimum influence of cold syringes on the temperature in the cell.

Preparing the Vapor Pressure Osmometer K-7000 for Measurement

Front View of the Instrument

You can locate the functional elements of the front panel with Fig. 1


Fig. 1 Front View of the Vapor Pressure Osmometer K-7000

Electrical Connections & Rear Panel

The vapor pressure osmometer can be used within a wide range of input voltages from 85 – 264 V, 47 – 440 Hz. The time-lag fuse is located below the socket for the power cable.

For problem-free operation connect the instrument only to sockets with a ground contact by using the three-pronged power cable included with delivery!

Rear Panel Elements


Fig. 2 Rear Panel of the Vapor Pressure Osmometer K-7000

Top View of the Instrument

The following functional elements are located on the top of the instrument, see Table 1 and Fig. 3.


Fig. 3 Top View of the Vapor Pressure Osmometer K-7000 (left) and Top View with Head Thermostat removed (right)

Table 1 Functional Elements of the Top of the instrument

| | |
|-----------------|--|
| Monitor Window | The window on the top of the osmometer allows to monitor the thermistors and to control the droplet sizes. The size of the droplets should be identical. The cell is illuminated by a small lamp and the light is deflected by a mirror to the window. |
| Syringe Ports | |
| Position Red | Pure solvent to rinse the right thermistor (reference) |
| Position Yellow | Pure solvent to rinse the left thermistor (sample) |
| Position 1 – 4 | For up to 4 different sample solutions |

Use Fig. 4 for identifying parts of the measuring cell.


Fig. 4 Elements and Design of Measuring Cell


Fig. 5 Side View with basic design of K-7000

Removing the Cell from the Instrument

Remove any syringes from the ports. Then lift the head of the osmometer. Unscrew the 2 female nuts, disconnect the black plug and lift the glass beaker including metal cover out of the thermostat. Remove the cover from the glass by slightly pulling the latter.

Cleaning the Thermistor Probes

Please rinse the thermistor probes which are located in the center of the cover with acetone followed by a solution of a surface active reagent and water. This procedure helps to adhere droplets on the thermistors.

Insertion of Vapor Wicks

Solvent wicks (order number A0429) are placed into the beaker to increase the liquid surface area for an optimum saturation of vapor in the entire cell. Proceed according to the following instructions: The large filter paper is rolled to a cylinder and is placed into the beaker with the window in the upper part. The small paper is rolled to a cylinder and the 2 paper rings are slipped on, one below and one above the window. The windows of both cylinders are adjusted to allow complete observation of the thermistors, see Fig. 6 next page.

Fill about 20 ml of the solvent for your measurements into the cell. The vapor wicks shall already be wetted during fill-up.

Insertion of the Measurement Cell

Place the beaker with the wicks into the thermostat. Both windows should face the front part of the osmometer. Put the cover onto the cell and reconnect the black plug.

Take care, that the plug is positioned properly! Otherwise the instrument will not work at all!

Be sure that all windows of the wicks are well aligned and face towards the front of the instrument, if not readjust. The ends of the 2 thermistor probes should be completely observable through the monitoring window. Reset the head thermostat with port # 1 in front.


Fig. 6 Preparing the vapor wicks

Operating your Vapor Pressure Osmometer K-7000

Switching the Vapor Pressure Osmometer K-7000 on

Use the on/off button on the rear panel side and switch on the instrument.

After power on the display shows the message **Knauer Osmometer** and a few seconds later the release number of the internal software. Startup is completed and the main menu will be displayed. The preselected temperature will be activated. All parameters of the last measurement will automatically be restored. However, you can select new parameters if desired.

Select a temperature $\leq 130^{\circ}\text{C}$, press key AUTOZERO, wait until the baseline is constant, then press the key TEST. After about 15 seconds either **Test ok** or **maintenance required** is displayed. In the last case the thermistor probe shall be sent back to the manufacturer for a routine check. It is recommended to run this test every 3 months or if measurements are not reproducible.

Using the Keys of the Vapor Pressure Osmometer K-7000

You can operate the instrument with the keys on its front panel, see Fig. 7. Table 2 gives a brief overview of the basic functions of these keys.


Fig. 7 Front Panel Elements of the Vapor Pressure Osmometer K-7000

Table 2 Basic Functionality of the Front Panel Keys

| | |
|----------------------|--|
| START STOP | To start and stop the measurement |
| AUTO ZERO | This function compensates the offset, finds the zero position if both thermistors hold droplets with pure solvent. You can use Autozero after 2 pure solvent droplets are attached and the baseline is constant. The measurement value will be reset to zero within 20 – 30 seconds |
| TEST | Automatic test run to check the performance of the thermistors |
| Right Arrow Key ➤ | Moves cursor to fields within a menu |
| Up Arrow Key ▲ | Selection of menu Increase of numeric values |
| Down Arrow Key ▼ | Selection of menu Decrease of numeric values |

Displays of the Vapor Pressure Osmometer K-7000

Main Menu

The main menu displays status information and allows access to other menus.

| | |
|---------------------|-------------------|
| SIGNAL:+0000 | 060°C +2°C |
| ◆ TIME:0.00 | Gain:064 |

Fig. 8 Example for a Main Menu Display

Table 3 Properties of the Display of the K-7000

| | |
|--|---|
| Signal/Measurement Value | +/- 4 digits |
| Page Selection | Black rhombus |
| Time | Measurement Time up to 30 min Seconds in units of 1/10 min |
| Selectable Gain ¹⁾ Settings | 1, 2, 4, 8, 16, 32, 64, 128, 256 |
| Selectable Cell Temperature | 20 – 130 °C |
| Selectable Δt ²⁾ | max. 6 °C |

¹⁾ The Gain Setting is related to the internal amplification of the raw signal and effects the output values of the instrument

²⁾ For the Vapor Pressure Osmometer K-7000 Δt is calculated as the difference between the Temperatures of Head and Cell $\Delta t = T_{\text{Head}} - T_{\text{Cell}}$

Selecting an appropriate Measurement Temperature

Prior to measurement you need to select the temperatures of the entire cell and the head individually. You can select both temperatures manually by selecting the appropriate field and edit the numeric values as follows. The head temperature is always entered as an increment of the cell temperature. If you are working with cell temperatures $\leq 35^\circ\text{C}$, set the head temperature zero.

1. Select Right Arrow Key ➤ until the corresponding position is reached.
2. Press either Up Arrow Key ▲ or Down Arrow Key ▼ to select the desired temperatures.

If you for example want to activate a temperature cell of 60°C and a head temperature of 62°C you can enter the values of 60°C and +2°C. The display shows a blinking value of the cell temperature until the preset temperature is actually reached.

To ensure stable temperature conditions, the cell temperature should be at least 5° above the room temperature.

Using the Time Display

Measurements can be started when the time display shows zero. The measurement time count can be started and stopped by pressing the start/stop key, respectively. This time window will automatically be stored unless not manually changed. You can use this feature for tracking the measurement time or certain intervals during measurements e. g. the signal reaches a consistent value. Pressing the start/stop button a second time will lead to a decremental time count until zero. On reaching zero the value that has been measured is kept in the display (HOLD position). Pressing the Start/Stop key again will lead to a decremental time count until zero. If the measurement is stopped during any time count, the value for time will automatically be reset to its initial value (no HOLD position!).

Identifying the Serial Number

Move the cursor using Right Arrow Key ➤ to the black rhombus and press Up Arrow Key ▲. In this screen the serial number of this particular instrument is displayed.

Fig. 9 Identifying the Serial Number on the display

Using the Calibration Mode

Press Up Arrow Key ▲ again. You should see a display similar to Fig. 10.


Fig. 10 Display using the mode calibration

You can use this feature for direct determination of osmolalities of physiological fluids.

For calibration purposes it is possible to change the measurement value to a lower number by reduction of the bridge voltage. We recommend values of $\geq 50\%$ compared to the initial value.

In the osmolality mode the measurement value can directly be correlated to an osmolality given in units of mOsmol/kg.

Example: Using the calibration mode

A standard solution of 400 Osmol/kg is measured. The gain is set to 16 and the measurement value is determined to be 455. In the calibration screen the actual value is 455.

Press the Right Arrow Key \blacktriangleright once. The cursor moves to the second line **set**.

Press the Down Arrow Key \blacktriangledown until 400 is reached. Press Right Arrow Key \blacktriangleright again: The calibration is now executed. In the display the word **calibrating** is shown. After about 10 seconds the actual value switches to 400. This change is due to an internal reduction of the bridge voltage. You can check the new voltage in the corresponding screen, see Fig. 11.

The calibration value is permanently stored.

The voltage applied to the Wheatstone Bridge is displayed in percent. Use this display to adjust and reset this value to 100%. During normal operation this value should give 100% unless a calibration has been performed and stored.

Using the Mode Actual Temperature

Use the Up Arrow Key \blacktriangle and select the following display. Use this screen for displaying the actual temperatures of cell and head.


Fig. 11 Temperature Display of Measurement Cell and Cell Head

Resetting the Bridge Voltage

To reset the bridge voltage to 100% press Right Arrow Key \blacktriangleright until the cursor reaches the field bridge value. Press Up Arrow Key \blacktriangle once. The display shows **Reset bridge voltage? Confirm by Cursor**. Press the Up Arrow Key \blacktriangle again and 100% are shown in the display.

Press Right Arrow Key \blacktriangleright first and Up Arrow Key \blacktriangle to activate the main menu.

Practical Choice of Solvents

Select the best suited solvent with the following criteria:

1. Good solubility of the sample.
2. No chemical interactions between sample and solvent.
3. Vapor pressure of sample is one tenth of vapor pressure of solvent or less.

Please find below influences of the solvent for your measurement:

4. The actual vapor pressure depression is related to the solvent used.
5. Especially with a large molecular mass, the sensitivity of the measurement is reduced because the number of moles in a solution compared to the mass becomes very small. Either high concentrated solutions have to be measured or very low particle concentrations must be accepted.
6. The solvent requires a good thermal stability in the selected temperature range.
7. The temperature difference between the measuring temperature and the boiling point of the solvent must be at least 15 °C.

Recommended measuring temperatures are listed in Table 4 "Some properties of most frequently used solvents for vapor pressure osmometry" on page 13.

Table 4 *Some properties of most frequently used solvents for vapor pressure osmometry*

| Solvent | Boiling Point °C | Measurement Range °C | Minimum Concentration mol/kg |
|-------------------------|---------------------|-------------------------|---------------------------------|
| Acetone | 56.2 | 37 - 45 | 0.005 |
| Acetonitrile | 81.6 | 37 - 60 | 0.010 |
| Benzene | 80.2 | 25 - 60 | 0.001 |
| Bromobenzene | 156.2 | 60 - 130 | 0.005 |
| Bromoform | 149.5 | 37 - 90 | 0.003 |
| Butyl Acetate | 126.1 | 45 - 90 | 0.005 |
| Carbon Disulfide | 46.3 | 25 - 37 | 0.001 |
| Carbon Tetrachloride | 76.6 | 25 - 45 | 0.001 |
| Chlorobenzene | 131.7 | 60 - 90 | 0.005 |
| Chloroform | 61.2 | 25 - 45 | 0.003 |
| Cyclohexane | 80.8 | 37 - 60 | 0.005 |
| Dichlorobenzene | 179.5 | 90 - 130 | 0.005 |
| Diethyl Ether | 34.6*) | 25 - 30 | 0.005 |
| N,N-Dimethyl Formamide | 155 | 90 - 120 | 0.005 |
| Dimethyl Sulfoxide | 100(d) | 60 - 90 | 0.005 |
| Dioxane | 101 | 60 | 0.005 |
| Ethanol | 78.3 | 45 - 60 | 0.002 |
| Ethyl Acetate | 77.1 | 37 - 60 | 0.002 |
| Ethyl Bromide | 38.3 | 25 | 0.005 |
| Ethylene Bromide | 131.6 | 60 - 90 | 0.005 |
| Ethylene Chloride | 83.6 | 25 - 60 | 0.003 |
| Formic Acid | 100.5 | 37 - 60 | 0.005 |
| Furane | 32*) | 25 | 0.001 |
| n-Heptane | 98.3 | 37 - 60 | 0.005 |
| n-Hexane | 68.8 | 37 - 60 | 0.005 |
| Mesitylene | 164.6 | 90 - 130 | 0.005 |
| Methanol | 64.7 | 45 | 0.005 |
| Methyl Bromide | 96.5 | 37 - 60 | 0.005 |
| Methyl Ethyl Ketone | 79.6 | 37 - 60 | 0.005 |
| Methyl-iso-Butyl Ketone | 116.9 | 45 - 60 | 0.005 |
| Methylen Chloride | 40.6*) | 25 | 0.003 |
| n-Propanol | 97.2 | 45 - 60 | 0.002 |
| Pyridine | 115.5 | 60 - 90 | 0.001 |
| Tertiary Butanol | 82.6 | 45 - 60 | 0.010 |
| Tetrahydrofurane | 65.5 | 45 | 0.005 |
| Toluene | 110.8 | 37 - 90 | 0.001 |
| Trichloro Ethylene | 87 | 37 - 60 | 0.001 |
| Water | 100 | 37 - 60 | 0.005 |
| o-Xylene | 143 | 90 | 0.001 |

- *) Low boiling solvents may evaporate. We recommend careful establishment of experimental conditions.
 (d) decomposition

Specifications of the required Syringes

Always use syringes which meet all specifications for the K-7000, see Fig. 12 (order number A0433)


Fig. 12 *Appropriate syringe for Vapor Pressure Osmometer K-7000*

This type of syringe consists of a glass cylinder, a metal piston with teflon head and a metal needle with a spring.

Performing Measurements

Using the syringes

Fill 2 syringes with pure solvent and place one in the red and the other in the yellow port. It is important to remove any air bubbles from the syringes prior to measurement. If the piston slips down by gravity it is recommended to empty the syringe and to press the teflon head on a plain surface to improve its sealing properties.

You can prepare 4 syringes with sample solutions. To avoid longer times for thermal equilibration, place all syringes in the desired ports as soon as possible. The syringe springs make sure that a defined distance from the thermistors is maintained. It is important that all ports are equipped with syringes, no matter if required for measurement.

To attach a droplet on a thermistor, push the syringe slightly down and thus move the needle towards the thermistor. By slightly pushing the piston downwards, droplets can be attached to the t probes. It is important to produce droplets with identical sizes in this step. Always rinse the thermistor 3 to 5 times with every sample before a final droplet is attached.

Finding the Zero Point

Prior to start a measurement you have to define the zero point of the instrument. Therefore both thermistors have to hold droplets of pure solvent.

1. Attach a droplet of pure solvent to both thermistor probes.
2. Wait for 5 minutes and make sure that a stable baseline is established.
3. Press button AUTOZERO.

Dropping Techniques

You can either exchange the sample droplet only or exchange sample and solvent droplet according to sections „Measuring new solvent and new sample droplets“ below and “Exchanging the sample droplet only” on page 15.

You can choose the desired technique according to your needs. It is important to decide for one of the techniques and constantly apply it. Measurement of new solvent and new sample droplets results in shorter measurement times.

Measuring new solvent and new sample droplets

To determine samples with different concentrations it is recommended to measure the samples in an order of increasing concentrations. To exchange both droplets in every measurement, proceed as follows:

1. Apply the sample droplet — The signal indicates a positive overload, see 1 in Fig. 13.
2. Apply the solvent droplet — The signal indicates a negative overload, see 2 in Fig. 13.
3. A constant measurement value is achieved when the thermal equilibrium is reached and the signal stops to increase, see 3 in Fig. 13.

The pattern of the signal may resemble graphs given in Fig. 13.


Fig. 13 Schematic Osmogram for a two droplets exchange

Exchanging the sample droplet only

1. Apply a new sample droplet — The signal indicates a positive overload, see 1 in Fig. 14.
2. A constant measurement value will be reached when the signal decreases until thermal equilibrium is reached, see 2 in Fig. 14.


Fig. 14 Schematic Osmogram for a one droplet exchange

Both schematic osmograms shown refer to head temperature setting of + 2°C with respect to the cell according to section "Selecting an appropriate Measurement Temperature" on page 10. With other temperature differences osmograms may differ from the given examples.

Thermistor Adjustment

If you cannot manage to attach droplets to any thermistor probe, proceed as follows.

1. Open the instrument.
2. Readjust the thermistors with careful manual bending until a distance of 0.9 cm is reached, see Fig. 4 on page 7. It is also possible to slightly bend the needle pipes.

Required Stabilization Time

Stabilization time is defined as the time interval required for the vapor pressure stabilization at the selected cell and head temperatures. Stabilization time begins when the preset temperatures are roughly reached and the respective values in the display stop flashing. Table 5 gives an overview of required stabilization times.

Table 5 Stabilization Times as a function of the measurement temperature

| Temperature [°C] | Stabilization Time [min] |
|---------------------|-----------------------------|
| 25 | 30 |
| 40 | 30 |
| 60 | 60 |
| 100 | 90 |
| 130 | 120 |

Three Different Measurement Modes

You can operate the Vapor Pressure Osmometer K-7000 in three different measurement modes.

Determination of the Total Osmolality of Physiological Fluids

When determining the total osmolality of physiological fluids, the measurement values are always proportional to the number of osmotically active particles within the measured fluid.

For calibration the K-7000 one solution of a standard with a known osmolality has to be measured. The calibration information is represented by the slope of the regression curve defined by the origin of the graph and the value $\frac{\text{Measured Value}}{\text{Known Osmolality}} \cdot K_{\text{calib}}$ unit is kg/mOsmol.

$$K_{\text{calib}} = \frac{MV}{c}$$

The osmolality of an unknown sample can be calculated with :

$$c \left[\frac{\text{mOsmol}}{\text{kg}} \right] = \frac{MV}{K_{\text{calib}}}$$

A commonly used standard is an aqueous solution of NaCl with 400 mOsmol/kg. Sets of ampoules are available (Order number Y1241). Anyway, other osmolalities can also be used (see table).


Fig. 15 Calibration Graph

Table 6 Relation between osmolality, freezing point depression and molality

| Real Osmolality $\frac{\text{mOsmol}}{\text{kg}}$ | Freezing Point Depression ΔT in $^{\circ}\text{C}$ | Concentration $\frac{\text{g NaCl}}{\text{kg H}_2\text{O}}$ | Molality $\frac{\text{mmol (Na}^+ + \text{Cl}^-)}{\text{kgH}_2\text{O}}$ | $R^*)$ |
|--|---|--|---|--------|
| 100 | 0.1858 | 3.09 | 105.6 | 1.06 |
| 200 | 0.3716 | 6.24 | 213.3 | 1.07 |
| 300 | 0.5574 | 9.45 | 323.1 | 1.08 |
| 400 | 0.7432 | 12.69 | 433.8 | 1.08 |
| 500 | 0.9290 | 15.94 | 545.0 | 1.09 |
| 600 | 1.1148 | 19.18 | 655.7 | 1.09 |
| 700 | 1.3006 | 22.45 | 767.5 | 1.10 |
| 750 | 1.395 | 24.10 | 823.9 | 1.10 |
| 1200 | 2.23 | 38.76 | 1325.1 | 1.10 |
| 1800 | 3.35 | 58.01 | 1983.2 | 1.10 |
| 2500 | 4.65 | 79.97 | 2734.0 | 1.09 |
| 3000 | 5.58 | 95.40 | 3261.5 | 1.09 |

*) R is molality divided by osmolality and thus represents the relative deviation of osmolality from molality

Note that it is extremely important to distinguish between osmolality and molality. In very diluted solutions both concentrations are comparable. But with increasing number of molecules in a solution the formation of aggregates is likely. The vapor pressure osmometer can only measure the real number of free particles, osmotically active particles in a solution (osmolality not molality). If aggregates are formed the osmolality is decreased compared with the molality.

Example: Determination of the total osmolality of a sample

Prepare the instrument as described in sections "Operating your Vapor Pressure Osmometer K-7000" on page 9 and "Selecting an appropriate Measurement Temperature " on page 10. Enter the following values:

- Cell Temperature T_C : 60 $^{\circ}\text{C}$
- Head Temperature T_H : 2 $^{\circ}\text{C}$
- Measurement Time: 1.5 min
- Gain: 16

Remember, the temperature is reached after the digits stop flashing, however the equilibration of the vapor pressure needs more time (Table 5 on page 16).

Performing a Calibration with NaCl of 400 mOsmol/kg

1. Open one ampoule (order number Y1241) of NaCl solution 400 mOsmol/kg,
2. Fill up one syringe and transfer it into syringe port 1. It is recommended to dispose the remaining residue of the ampoule. Use it only immediately after the ampoule is opened.
3. Fill your sample solution into another syringe and transfer it into port 2.
4. Fill distilled water into two further syringes and place them in the yellow and red ports, respectively.
5. Insert empty syringes into the remaining ports 3 and 4 to maintain temperature consistency.
6. Make sure that prior to the first measurement the stabilization time, see Table 5 Stabilization Times as a function of the measurement temperature, on page 16, is considered.
7. Solvent drops are attached to either thermistors. After the baseline is on a constant level the button Autozero is pressed. The measurement value shown in the display will turn to zero after 20 – 30 seconds. The osmometer is ready for measurement.
8. Rinse the left hand thermistor with 3 – 5 droplets of the standard solution, the last drop should remain at the end. It is important to apply droplets of an equal size.
9. Press the button START immediately after the droplet is attached.

The measurement signal shows **OVL** which means that the very sensitive thermistor realizes a temperature change. After 20 – 30 seconds the signal is back in the range and approaches a constant value. After 1.5 min the display shows the measured value in the **HOLD** position.

This standard has to be measured at least 3 times. The standard deviation should not exceed 0.5%.

Example: Evaluation of calibration data

You may find data from your calibration that can be evaluated according the following example:

Table 7 Example of Calibration Data

| Number of injection | c in $\frac{\text{mOsmol}}{\text{Kg}}$ | Measurement Value |
|-----------------------------|--|-----------------------------------|
| 1 | 400 | 43.37 |
| 2 | 400 | 43.51 |
| 3 | 400 | 43.47 |
| 4 | 400 | 43.61 |
| 5 | 400 | 43.54 |
| Mean | | 43.50 |
| Standard Deviation | | 0.09% |
| $K_{\text{calib}} = 0.1088$ | | $\frac{\text{kg}}{\text{mOsmol}}$ |

Evaluating Measurement Data

Rinse the left hand thermistor with 3 – 5 droplets of distilled water. The last droplet should remain at the end of the thermistor. Check the zero level.

Now measure your sample at least 3 times as described above.

Table 8 Example of Measurement Data and Evaluation based on calibration

| Number of injection | Measurement value |
|---|--|
| 1 | 34.43 |
| 2 | 34.61 |
| 3 | 34.50 |
| | Mean 34.51 |
| | Standard Deviation 0.03% |
| Total osmolality of the sample : | |
| | $34.51 \frac{\text{mOsmol}^*)}{0.1088 \text{ Kg}} = 317 \frac{\text{kg}}{\text{mOsmol}}$ |

*) Taken from calibration, see Table 7

Molecular Mass Determination for M < 500 g/mol

The number of particles or moles of a substance dissolved in a solution influences the measurement value in osmometry. The molecular mass can be determined if the sample concentrations are known. This is only valid for diluted solutions.

The basic relation for molecular mass determination equals:

$$\frac{n \text{ Sample}}{m \text{ Solvent}} = \frac{m \text{ Sample}}{m \text{ Solvent} * M \text{ Sample}}$$

with

- n = Number of sample molecules
- m = Mass of sample and solvent respectively
- M = Molecular mass of the sample:

For a molecular mass less than 500 g/mol the measurement value is proportional to the number of moles. This is very similar to the osmolality determination.

A sample is measured where concentration and molecular mass are known (c in mol/kg). The slope of the graph 0/0 and measurement value/concentration is K_{calib} with kg/mol.

$$K_{\text{calib}} = \frac{MV}{c}$$

The unknown molecular mass of a sample can be determined with a known concentration of the unknown substance given as c with g/kg. The slope of the regression curve passing the origin of the graph and the measured value K_{meas} with kg/g. The molecular mass of the unknown sample is determined with:

$$M \left[\frac{\text{g}}{\text{mol}} \right] = \frac{K_{\text{calib}}}{K_{\text{Meas}}}$$


Fig. 16 Calibration Graph

Example: Calibration and Determination of Molecular Mass

The calibration is performed out with benzil. An example for calibration data is given in Table 9. Table 10 contains example data for real measurements and for the calculation of the molecular mass using calibration and measurement data.

Table 9 Calibration with Benzil

| Number of injection | c in $\frac{\text{mol}}{\text{Kg}}$ | Measurement value |
|---------------------|-------------------------------------|---|
| 1 | 0.042 | 22.35 |
| 2 | 0.042 | 22.19 |
| 3 | 0.042 | 22.24 |
| Mean | | 22.26 |
| Standard Deviation | | 0.04% |
| | | $K_{\text{calib}} = 530 \frac{\text{kg}}{\text{mol}}$ |

Table 10 Sample Measurement

| Number of injection | c in $\frac{\text{g}}{\text{Kg}}$ | Measurement value |
|------------------------|-----------------------------------|---|
| 1 | 18.45 | 31.41 |
| 2 | 18.45 | 31.52 |
| 3 | 18.45 | 31.45 |
| Mean | | 31.46 |
| Standard Deviation | | 0.02% |
| | | $K_{\text{Meas}} = 1.71 \frac{\text{kg}}{\text{g}}$ |
| Molecular Mass: | | $M = \frac{530}{1.71} = 309.9 \text{ g/mol}$ |

Molecular Mass Determination for $M > 500$ g/mol

For larger molecules the relation between measurement value and concentration is not linear anymore. Therefore, the molecular mass has to be determined as follows:

Calibration

1. Select a standard with a molecular mass in the range of the expected ones of your samples.
2. Prepare 3 – 5 different concentrations
3. Measure these standards.
4. Determine the quotients of measured values and the corresponding concentrations and plot these data vs concentration (c with mol/kg).
5. Determine the regression curve and extrapolate to concentration $c = 0$ for K_{calib} in, see below...


Fig. 17 Calibration Graph for $M \geq 500$ g/mol — Example Data

Measurement:

1. Prepare 3-5 concentrations of each sample
2. Measure these samples.
3. Plot measurement value/concentration vs concentration c .
4. Find the regression curve and extrapolate to concentration $c = 0$ for K_{meas} determination in kg/g is determined, see below.


Fig. 18 Measurement Graph for $M > 500$ g/mol — Example Data

The molecular mass can be determined with

$$M \left[\frac{\text{g}}{\text{mol}} \right] = \frac{K_{\text{calib}}}{K_{\text{Meas}}}$$

Determination of the average molecular mass of a Polystyrene

Table 11 Calibration with PS 10100

| Number of injection | c in $\frac{\text{mol}}{\text{kg}}$ | Measurement Value MV | $\frac{MV}{c}$ in $\frac{\text{kg}}{\text{mol}}$ |
|---------------------|-------------------------------------|----------------------|--|
| 1 | 0.001 | 0.593 | 593 |
| 2 | 0.002 | 1,28 | 640 |
| 3 | 0.003 | 2.13 | 710 |

Extrapolation to c = 0 gives: $K_{\text{calib}} = 527.6 \text{ kg/mol}$

Table 12 Measurement of an unknown polystyrene sample

| Number of injection | c in $\frac{\text{g}}{\text{kg}}$ | Measurement Value MV | $\frac{MV}{c}$ in $\frac{\text{kg}}{\text{g}}$ |
|---------------------|-----------------------------------|----------------------|--|
| 1 | 11.50 | 0.71 | 0.0617 |
| 2 | 20.12 | 1.25 | 0.0621 |
| 3 | 28.79 | 2.07 | 0.0719 |

Extrapolation to c = 0 gives: $K_{\text{Meas}} = 0.054 \text{ kg/g}$

Molecular mass of the sample:

$$M = \frac{527.6}{0.054} = 9770 \text{ g/mol}$$

Operating the K-7000 with EuroOsmo® 7000 Software

The full capabilities of the K-7000 vapor pressure osmometer are accessed under operation with the EuroOsmo® 7000 software for Windows, version 1.2 or higher.

The License Agreement on page 69 should be read before using the EuroOsmo®7000 software! Using the software, the user indicates that he has read the License Agreement and agrees with the regulations.

This chapter will explain the basic knowledge needed for working with K-7000 and the EuroOsmo® 7000 software.

Hardware Requirements and Installation of the EuroOsmo® 7000 Software

| | | |
|-------------------------------|---|-------------------------|
| CPU | Pentium 130 MHz | |
| Memory (Minimum requirements) | | |
| | Windows NT 4.0 | Windows 95 |
| O/S | 32 MB Service Pack 3 (or later) (The Service Packs must have the same language as the O/S!) | 16 MB Service Pack 1 |
| Disk | 20 MB free disk space | |

Before installation make sure you have the appropriate version of the Operating System. If not, you must install the corresponding Service Pack.

Installation of the EuroOsmo® 7000 software from a CD ROM

1. Insert CD into your CD drive
2. Wait up to 30 seconds. The set-up starts automatically
3. If your CD ROM drive does not have AutoPlay:
 - Double-click on „My Computer“
 - Right-mouse click on your CD ROM
 - Select AutoPlay from the right-click menu
 - If there is no AutoPlay menu, select Explore from the right-click menu, go to „EuroOsmo 7000/Disk1“ and double-click on „Setup.exe“
4. Follow the installation guide.

Installation of the EuroOsmo® 7000 software from a floppy disk

1. Insert the Disk #1 into your drive
2. Double-click on „My Computer“
3. Double-click on „3.5 Floppy (A:)“
4. Double-click on „Setup.exe“
5. Follow the installation guide.

Hardware installation

The K-7000 vapor pressure osmometer has to be connected with a defined port of the computer (com 2 for instance) prior to starting the program.

Explanation of the main screen

After starting the program by a double-click on the corresponding icon or with START / PROGRAMM / EUROOSMO you will get a screen similar to Fig. 19.


Fig. 19 Opening screen of the EuroOsmo® 7000 software

Table 13 Description Screen Elements

| Menu line | File Edit View Run Help |
|--------------|---|
| File | Menu for opening, saving, and printing of data files or previews (measurement or calibration) |
| Edit | Menu for editing file changes and search performing |
| View | Menu to show or to hide status bar and toolbar |
| Run | Remote start and fixed gain can be set and Autozero can be performed |
| Help | The Version of the EuroOsmo® is shown |
| Toolbar | |
| | Most of the icons are usual windows features and do not need further explanation |
| | |
| | Paste the data of a selected (first) field into following lines of the chosen column |
| | |
| | perform an Autozero |
| Working mode | Select the type of measurement |
| Osmol | Determination of the total osmolality of physiological fluids |
| MW < 500 | Determination of molecular mass up to 500 g/mol |
| MW > 500 | Determination of molecular mass exceeding 500 g/mol |


Fig. 20 Upper monitor part (measurement view) depending on the operation mode

Table 14 Description Status Field

| | |
|-------------------------|--|
| System Status: | The data are taken from the setup (below). The temperature values are blinking, as long as the actual values differ from the settings. |
| Sample Number: | Selected measurement series according to the Sample List Table (see below). All the graphic windows will be updated corresponding to the selected sample. A double click on a sample list row is in the same manner effective. |
| Start: | Button to start a measurement or calibration run |
| t: | time after which the measurement value is acquired and the results are calculated |
| Scroll/ ... / ... / ... | Refers to the monitoring part on the right-hand side in either the measurement or calibration view. Use the < and > buttons to scroll. |

You have the choice between four different view options in the upper right part of the screen:

1. RT Monitor (real time monitor)
2. Measurement View
3. Calibration View
4. Full View

The RT monitor serves to check the base line. The gain here can be changed by using the or buttons. The RT monitor is self-scaling. The scaling is based on the maximum signal within the actual time window.

The measurement view and the calibration view do not differ from each other. Both of them show four of up to 20 possible runs of a sample monitored.


Fig. 21 Monitoring of calibration runs

Please use the < and > buttons to scroll the runs. The control fields below each of the curves can be used to accept or reject single runs. The gain and signal values and as far as for the measurement mode necessary the concentrations are also given and scrolled in this view (compare also Fig. 20 on page 25). Only in case of determination of molecular masses below 500 g/mol (B), the measurement and calibration views differ from each other. Its calibration view looks like the measurement view of osmolality mode (A). The time of the measurement value acquisition is marked by a red triangle. Any measurement or calibration curve will be enlarged in the full view up to this triangle by a double click on the corresponding curve window. In the upper right corner you it is indicated, which curve is selected. It is the calibration curve no. 2 in the shown example marked by a C. A leading letter M indicates measurement curves.


Fig. 22 Full view example of a calibration curve

For any of these four views the graph settings can be changed, if the Setup is selected in the lower part of the screen. On the right hand side you can select appropriate colors for the background, the axis, and the curves as well. Also you can hide or show the grid by clicking the control field, see Fig. 23.

Software set up for the measurement


Fig. 23 Setup screen

First of all you have to determine the communication port for the connection to your K-7000. This is only necessary when you start your osmometer K-7000 with EuroOsmo® for the first time. If you intend to change the port for any reason, this change becomes effective only by restarting the EuroOsmo® software. Then you have to select the working temperature (20-130°C), the temperature difference between cell and head thermostat (0-6°C), the measurement time, and the RT monitoring time.

The measurement time sets the time window for plotting the graph. It is not identical with timer value *t*, determining the time period after which the measurement value is acquired. The measurement time has to be larger than this timer value. Consequently the RT monitoring time should be larger than the measurement time.

Further settings are the used solvent, the operator's name, and last not least the start gain. During the measurement the software selects automatically the optimum gain depending on the signal height. Otherwise a fixed gain can be set in the RUN menu. The gain is relevant for the curve plots. The measurements are performed with the maximum gain.

All the settings in this field will become effective to the selected measurement mode by clicking the **Apply** button.

Both the other possibilities of the lower screen part (**Calibration** and **Sample List**) will be explained in detail with the description of the measurement performance in the three different modes.

Performing measurements

For the following sections please compare with the corresponding chapters where the measurements without EuroOsmo® software are described. Here only the changes are described, due to the use of the software package. Common characteristics of K-7000 as for instance the required warm up time or the technique of droplet exchange will not be mentioned here. Before starting **any** measurement or calibration, you have to determine the zero point by performing an AUTOZERO. This is independent on the kind of measurement. For this purpose proceed as described in the section "Finding the Zero Point" on page 14 and click on the  button instead of pressing the AUTOZERO button at the K-7000.

Determination of the total osmolality

Calibration

Choose the working mode Osmol and carry out the set up as described above. Select the calibration view in the upper right part of the screen as well as calibration in the bottom part. Enter the name of the calibration substance and the osmolality of your standard solution.

After exchanging the solvent droplet on the sample thermistor (left) by a solution droplet click immediately on the **Start Calibration** button on the right hand side of the lower calibration window. The measurement will be automatically stopped when the selected time **t** has passed. The applied gain and the measured value are listed below the osmogram (e.g. 5: 4/1.1).

The value of K_{calib} is calculated and shown. It is recommended to carry out at least 3 measurements to receive reliable values. In case that any of the calibration runs will fail you can omit this run by erasing the accept tick 4 in the corresponding control field.

Measurement

Switch to the Measurement View and Sample List. Sample number "1" is automatically selected. To confirm the calibration it is recommended to measure the NaCl standard solution as a sample again. After the droplet is attached to the sample thermistor (left) please start the measurement immediately by clicking the **Start** button in the central window on the left side of the screen. The measurement will be finished after the selected time period has expired. The osmolality of the sample will be automatically calculated and listed in the sample table.

You can measure each sample up to 20 times. All values of this sample are calculated to one mean value of osmolality. Again, if any run has failed – e.g. the droplet may be fallen down from the thermistor – you can omit this value simply by erasing the accept tick 4 in the corresponding control field. The value in the sample list represents the mean value of all accepted runs. The #-field indicates the quantity of the accepted runs.

Determination of molecular mass <500

Select the working mode **MW<500** and perform the set up and the AUTOZERO as described above.

Calibration

Select the calibration view in the upper right part of the screen as well as calibration in the bottom part.



Fig. 24

Calibration field for determination of molecular masses <500

Enter the calibration substance and the concentration of your standard. The concentration can be entered directly in mol/kg, or by clicking on the button (arrow marked in Fig. 24) an additional window for concentration calculation will emerge. With this second option the mass of the standard, its molecular

mass and the mass of the solvent can be entered and the concentration will be calculated automatically.

The image shows a dialog box titled "Concentration Calculation" with a close button (X) in the top right corner. It contains four input fields: "Sample Weight (g)", "Sample Molecular Weight (g/mol)", "Solvent Weight (kg)", and "Sample Concentration (g/kg)". To the right of these fields are two buttons: "Apply" and "Cancel".

Fig. 25 Concentration calculation window

For further proceeding the calibration continue as in the osmolality mode. It is strongly recommended to carry out at least 3 measurements to receive reliable values. Furthermore we advise to perform a test run with the calibration solution as a sample.

Measurement

The measurement runs are carried out in the same way as for the determination of the osmolality. The only difference is to enter the concentration of your sample solution. Click on the button on top of the measurement view Fig. 26 to open the concentration calculation window similar to that in Fig. 25. In this window, there is no access to the molecular mass field.


Fig. 26 Measurement view for determination of molecular masses <500

The measurement results will be monitored again in the sample list. The shown values represent the mean molecular masses and the #-field indicates the number of accepted runs.

Determination of molecular mass >500

Select the measurement mode **MW>500** and perform the set up and the AUTOZERO as described above.

Calibration

Select the calibration view in the upper right part of the screen as well as calibration in the bottom part. Enter the calibration substance in the bottom part of the screen. The individual concentrations of your standard solutions must be entered for every single run above the corresponding osmogram windows. For this purpose open the concentration calculation windows for each run you are going to perform.

After finishing the second calibration run, a curve similar to Fig. 17 on page 21 will be shown in the calibration view at the bottom of the screen, see Fig 27.


Fig. 27 *Extrapolation curve of measured calibration values*

This curve will be completed with every new or additional calibration run. If you run a given concentration twice or more, the software automatically calculates the mean value of these runs and determines K_{calib} by extrapolation of the graph $MV/c = f(c)$ towards $c \rightarrow 0$.

Measurement

For the measurement proceed in a similar way. Again the concentration must be entered separately for all runs. As soon as at least two runs are carried out a curve similar to Fig. 18 on page 21 will emerge in the calibration view at the bottom of the screen. You can swap between the calibration graph (Fig. 27) and the measurement graph by clicking on the measurement control field.


Fig. 28 *Extrapolation curve of measured values*

For each sample from the sample list you can perform up to 20 runs at up to five different concentrations. The software automatically calculates the mean values for each measured concentration and extrapolates the K_{meas} similar to K_{calib} . The calculated molecular mass is also indicated in the sample list.

Data saving

All data achieved by the EuroOsmo® software can be saved in sample or method files. Any of these files include the data of

1. method set up
2. calibration
3. measurement
4. evaluation

The files are characterized by the extension **.osm**. The files will be opened and closed in the convenient manner by using the window's icons of the toolbar.

At any time only one of these files can be opened. Therefore, if you are going to open a given file you always will be asked for saving the already opened (untitled) one. Consequently you also will be asked for saving the changes in file when you close the EuroOsmo® program.

For a multiple use of calibration data these can be saved in separate calibration files with the extensions

- .cls** osmolality;
- .cli** molecular masses < 500
- .clg** molecular masses > 500.

These files are opened and saved by using the respective buttons on the bottom part of the screen after selecting the calibration table.


Fig. 29 Opening and saving of calibration files

Trouble Shooting Tips

| Problem | Probable Cause and Advice |
|---|---|
| The instrument is not working | Check the correct positioning of the electrical connector at the cell cover. |
| The Instrument remains after conditioning in overflow | The set temperature is too low. It should be at least 5° above room temperature. |
| Needles of the syringes always pass the thermistors | Thermistors need to be readjusted. The distance between both probes should be 0.9 cm. You can also adjust the directional pipes for the needles. |
| Drops always fall down and cannot be permanently attached to the thermistor | Rinse the thermistors with acetone followed by a surface active reagent and water. |
| Zero line not adjustable | Defective thermistors, please press the test button (proceed as in Operating your Vapor Pressure Osmometer K-7000 on page 9) |
| Baseline drift | Temperature equilibrium not yet reached, wait some time before continuing |
| Syringes bleed | Remove the piston and press it with the teflon end on a table to reset the seal to the glass cylinder |
| Measurement values are not reproducible | Temperature equilibrium not yet reached, wait some time before continuing |
| No constant measurement value | Droplet size is not consistent, the smaller the droplet the larger the measurement value. |
| | Equilibration time too long, use a higher temperature or a different solvent |
| Measured Value: Increase followed by an immediate decrease | Volatile sample, cannot be measured |
| Measurement graphs have different shape | Decomposition of the solvent, please exchange |
| Measurement graphs have different shape | Dropping technique is not consistent |

Spare Parts & Accessories

Table 15 List of Standard Supplies (included in Package A3701)

| | |
|----------|--|
| A3703 | Universal Thermistor Probe (20 – 130 °C) |
| A0807 | Glass Beaker with Glass Cylinder |
| A0429 | Package with 10 Vapor Wicks |
| 6x A0433 | Syringe with Teflon/Steel Piston |
| Y1277 | Cleaning Solution |
| Y1026 | 10 g Benzil |
| | Recorder Cable |
| | Power Cable (110 V or 220 V) |
| V7109 | Operation Manual |

Table 16 List of additional accessories

| | |
|-------|--|
| A0428 | Syringe with Glass Piston |
| A0430 | Package with 12 Needles and 12 springs |
| Y1241 | Package with 10 Ampoules NaCl Calibration Solution c=400 mOsmol/kg |
| A3702 | EuroOsmo [®] 7000 Software Package for Windows 95 / 98 / NT, (RS232 cable included) |

Data Sheet of the Vapor Pressure Osmometer K-7000

| | |
|--------------------------|---|
| Features | Determination of the average molecular mass of polymers in aqueous and organic solutions Determination of the total osmolality of physiological fluids External Data Evaluation with EuroOsmo 7000 Software Package |
| Molecular Mass Range | Aqueous: $\leq 10,000$ g/mol Organic: 40 – 40,000 g/mol |
| Concentration range | 0.001 – 15 molal |
| Lower Detection Limit | 3.3×10^{-5} mol/kg in toluene 1.7×10^{-4} mol/kg in water |
| Temperature Consistency | $\pm 1 \times 10^{-3}$ °C |
| Cell Temperature Range | 20 – 130 °C |
| □T Head Thermostat | max. 6 °C |
| Warm Up Times | □ 40 °C: 0.5 h □ 60 °C: 1.0 h □ 100 °C: 1.5 h □ 130 °C: 2.0 h |
| Syringes | 1 ml. glass cylinder with teflon-headed steel piston |
| Min. Cell Solvent Volume | 20 ml |
| Gain | 1 – 256 |
| Voltage Range | 85 – 264 V / 47 – 440 Hz. 60 VA |
| PC Connection | RS 232 |
| Recorder Output | 1 V |
| Dimensions | 160 x 182 x 340 mm |
| Weight | 5.4 kg |

We reserve the right to make technical changes for instrument improvement. Minor changes will be self-explanatory. Improvements affecting use or maintenance will be described in supplementary pages to this manual.

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Warranty statement

The Vapor Pressure Osmometer K-7000 is guaranteed to be of quality material and workmanship. This warranty is valid for one year from the date of shipment. Warranty is void by improper operation, application and/or maintenance of the instrument. Any expendable components such as glass parts and thermistors are excluded from this warranty. This warranty is validated when, within fourteen days of receipt, the accompanying guarantee card is completed and returned to the manufacturer:

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If we find a defect covered by the guarantee, repair or replacement, at our discretion, will be carried out free of charge. Packing and transport costs are borne by the purchaser.

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Declaration of conformity

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Manufacturer's name and address:**Herstellername und -adresse**

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Dampfdruckosmometer K-7000

Order Number, Bestellnummer **A3701**

complies with the following requirements and product specifications:

- Low Voltage Ordinance (73/23/EWG)
EN 61010-1 (08/2002)
- Engineering Guidelines (89/392/EWG)
- EMC Ordinance (89/336/EWG)
EN 6100-3-2 (10/2006)
EN 61326-1 (10/2006)

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KNAUER Software License Agreement

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(hereinafter Licensor or KNAUER)

and the software user (hereinafter Licensee)

for the software package

EuroOsmo[®]7000 V1.2 and handbooks.

In return for acquiring a license to use the software ("EuroOsmo[®]7000 V1.2") and related documentation, Licensee agrees to the following terms and conditions:

1. License.

This agreement grants the Licensee a license to

- (a) use the EuroOsmo[®]7000 V1.2 on a single computer system or, in the case of a multi-user or networked system which permits access to the EuroOsmo[®]7000 V1.2 by more than one user at the same time, at one (1) single working location, and
- (b) make one copy of the EuroOsmo[®]2000 V1.2 in machine readable form solely for back-up purposes provided your reproduce KNAUER's copyright proprietary legends and mark the back-up copy with the express wording "back-up copy Knauer GmbH, Hegauer Weg 38, D-14163 Berlin".

2. Restrictions.

Licensee may not distribute copies of the EuroOsmo[®]7000 V1.2 to others or electronically transfer the EuroOsmo[®]7000 V1.2 from one computer to another over a network. Licensee may not use the EuroOsmo[®]7000 V1.2 from multiple locations of a multi-user or networked system at any time. The EuroOsmo[®]7000 V1.2 contains trade secrets and in order to protect them the Licensee may not decompile, reverse engineer, disassemble, or otherwise reduce the Software to a human-perceivable form: LICENSEE MAY NOT MODIFY, ADAPT, TRANSLATE, RENT, LEASE, LOAN, RESELL FOR PROFIT, DISTRIBUTE, NETWORK OR ANY DERIVATIVE WORKS BASED UPON THE SOFTWARE OR ANY PART THEREOF. The Licensee must not remove from the program firm names, trademarks, logos, copyright notices and other markings or devices, which serve the protection of the program.

3. Ownership of EuroOsmo[®]7000 V1.2.

Licensee owns the media upon which the EuroOsmo[®]7000 V1.2 is recorded or fixed, but KNAUER, Licensor, retains title and ownership of the EuroOsmo[®]7000 V1.2 recorded on the original media and all subsequent copies of the EuroOsmo[®]7000 V1.2, regardless of the form or media in which or on which the original and other copies may exist. This license is not a sale of the EuroOsmo[®]7000 V1.2 or any copy.

4. Confidentiality.

Licensee agrees to maintain the EuroOsmo[®]7000 V1.2 in confidence and to not disclose the EuroOsmo[®]7000 V1.2 to any third party without the express written consent of KNAUER. Licensee further agrees to take all reasonable precautions to preclude access of unauthorized persons to the EuroOsmo[®]7000 V1.2. The Licensee must not pass copies of the program or the handbook to third persons.

5. Term.

This license is effective until January 1, 2020, unless terminated earlier. Licensee may terminate the license at any time by destroying the Software (including the related documentation) together with all copies or

modifications in any form. KNAUER will have the right to terminate the Licensee's license immediately if the Licensee fails to comply with all copies or modifications in any form. The termination does not raise Licensee's claim for refunding of the license fee.

6. Limited Warranty.

KNAUER warrants only that the media upon which the EuroOsmo®7000 V1.2 is furnished will be free from defects in material or workmanship under normal use at the time of delivery. KNAUER DOES NOT AND CANNOT WARRANT THE PERFORMANCE OR RESULTS LICENSEE MAY OBTAIN BY USING THE SOFTWARE OR DOCUMENTATION. THE FOREGOING STATES THE SOLE AND EXCLUSIVE REMEDIES KNAUER WILL PROVIDE FOR BREACH OF WARRANTY. EXCEPT FOR THE FOREGOING LIMITED WARRANTY, KNAUER MAKES NO WARRANTIES; EXPRESSED OR IMPLIED, AS TO NON-INFRINGEMENT OF THIRD PARTY RIGHTS, MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. Some countries do not allow the exclusion of implied warranties or limitations on how long an implied warranty may last, so the above limitations may not apply to the Licensee. This warranty gives the Licensee specific legal rights and the Licensee may also have other rights which vary from country to country.

7. Limitation of Liability.

IN NO EVENT WILL KNAUER BE LIABLE TO LICENSEE FOR ANY SPECIAL DAMAGES, INCLUDING ANY LOST PROFITS, LOST SAVINGS, OR OTHER INCIDENTAL OR CONSEQUENTIAL DAMAGES - EVEN IF KNAUER HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES, OR FOR ANY CLAIM BY ANY OTHER PARTY. Some countries do not allow the exclusion of limitation of special, incidental, or consequential damages, so the above limitation, or exclusion may not apply to the Licensee. Any liability is – in any case – limited to the license fee.

8. Limitation of Remedies.

KNAUER's entire liability and Licensee's exclusive remedy shall be:

- (a) the replacement of any media not meeting KNAUER's limited warranty which is returned to KNAUER, or
- (b) if KNAUER or its distributor is unable to deliver replacement media which is free of defects in materials or workmanship, the Licensee may terminate this Agreement by returning the EuroOsmo®7000 V1.2 and the Licensee's money will be refunded.

9. This license is not assignable.

The Licensee does not have the right to assign rights or grant sublicenses.

10. License Fee / Royalty.

For the license to use the contractual EuroOsmo®7000 V1.2 in the above defined scope and within the restrictions the Licensee pays to the Licensor a fee. The amount becomes due and payable upon delivery of the program discs. The license to use of the program commences with payment of the royalty.

11. If a provision of this Agreement should be invalid or become invalid, then the legal effect of the other provisions shall remain unchanged. Instead of an invalid provision, a valid provision shall be deemed to be agreed which comes closest to what the parties intended economically.

This document is the complete Agreement between KNAUER and Licensee; it supersedes all and any prior agreement or offer. All and any modifications of this agreement have to be made in writing.

For all controversies out of this Agreement, the District Court Berlin shall have exclusive jurisdiction. The applicable law shall be that of the Federal Republic of Germany (German Civil Code) excluding the UN-Code on the international sale of goods.