Smartline Valve Drives S6, S12, S16 Manual / Handbuch V7447-1 09/2007

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CONTENTS

Using this Manual	4
Conventions in this manual.....	4
SOP's in this manual	4
General Description.....	5
Preparing the Electric Valve Drives S 6/-12/-16 for Operation.....	5
Standard Delivery List	5
Front view of the Valve Drives S 6/-12/-16.....	6
Rear view of the Valve Drives S 6/-12/-16	7
Operating the Valve Drive S 6/-12/-16	8
Power supply	8
Mounting the valve.....	8
Configuring the Smartline Electric Valve Drive	8
6-Port-Multi Channel Valve Switch Positions :	10
Controlling the Valve Drives S 6 /-12 /-16	11
Manual control	11
External Control of the Valve Drives S 6/-12/-16.....	11
Assembling a remote control cable.....	14
Valve drives S 6/12/16 as parts of the Smartline Assistant 6000.....	14
Troubleshooting.....	16
Spare parts and accessories	16
Maintenance by the user	16
Declaration of conformity.....	30
Guarantee statement.....	31

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Using this Manual

This manual refers to the Smartline Electric Valve Drives S 6, S 12, and S 16. It is valid for any combination of the drive S 6 with the 6-port/3-channel-injection valve, 6-port/2-channel-switching valve, 7-port/1-channel- switching valve, the drive S 12 with the 13-port/1-channel- switching valve, and the drive S 16 with the 17-port/1-channel- switching valve.

Conventions in this manual

Important Hints are marked by the marginal hand symbol.

Special warnings are indicated by the marginal warning sign and printed in bold letters.

The marginal lamp symbol indicates **helpful advice's**.

SOP's in this manual

The **Standard Operating Procedures (SOP)** provided with this manual offer a convenient way of structuring complex tasks for operating the Smartline Electric Valve Drive S 6, S 12, or S 16. They include step-by-step instructions assisting the user through all routine tasks during operation. They can be used for documentation purposes: They can be copied, applied, signed, and filed to document the performance of the instrument.

It is very important to follow all instructions and SOP's in this manual in order to operate the instrument and accessories. This ensures proper results and longevity of your equipment.

SOP 1	Mounting the valve to the electric valve drive	8
SOP 2	Connecting the plug strips	14

General Description

The Smartline Electric Valve Drives S 6, S 12, and S 16 are in design and operation for the most part identical. They are devised for automatic control of the 6-port/3-channel-injection valve, 6-port/2-channel-switching valve, 7-port/1-channel-switching valve by S 6 drive, the 13-port/1-channel-switching valve by S 12 drive, and the 17-port/1-channel-switching valve by S 16 drive respectively. The valve drive has to be adapted to the used valve by means of a set of DIP switches at the bottom of the drive housing.

All KNAUER valve drives can be integrated into the Smartline Assistant 6000.

Valves with 2 positions

The 6-port/3-channel-injection valve and the 6-port/2-channel-switching valve can be switched between the LOAD and INJECT position. Switching the positions can be done with a 60° turn forwards or backwards.

Valves with 6 positions

The 7-port/1-channel-switching valve can be switched to 6 positions. Each switch command will result in a 60° turn, so that the 6th switch command returns to the first position.

Valves with 12 positions

The 13-port/1-channel-switching valve can be switched to 12 positions. Each switch command will result in a 30° turn, so that the 12th switch command returns to the first position.

Valves with 16 positions

The 17-port/1-channel-switching valve can be switched to 16 positions. Each switch command will result in a 22.5° turn, so that the 16th switch command returns to the first position.

Using remote or software control valves with more than two positions can be switched to any one directly.

Preparing the Electric Valve Drives S 6/-12/-16 for Operation

Unpacking

After unpacking, please check the device and accessories thoroughly for any damage that may have occurred during transportation. If necessary, contact the carrier to claim any damages.

Use the "Standard delivery list" to check if the electric valve drive is complete. Please contact our service department if anything is missing or if you need support. It is important to fill out the guarantee card and return it to us as soon as possible.

Please remove the protection foil from the display!

Standard Delivery List

1	Smartline Electric Valve Drives S 6, S 12, or S 16
1	Power supply unit 24V/50W
1	Power supply cable, 230V
1	Adapting plate 70x50 Al
1 Set	Tools
1	RS-232 cable
1	WAGO Terminal Strip
1,5 m	Set signal conductors
1	User Manual

Front view of the Valve Drives S 6/-12/-16


Fig. 1 Front view of the Electric Valve Drive S 16 with 17port/1channel-switching valve

- 1 Display
- 2 Foil keys
- 2.1 LOAD position, or increasing the position number
- 2.2 INJECT position, or decreasing the position number
- 3 Valve (7port/1channel)
- 4 Adapter plate
- 5 Mounting screw


Fig. 2 Schematic front view of the 7-port/1-channel-switching valve


Fig. 3 Schematic front view of the 6-port/2-channel-switching valve and the 6-port/3-channel-injection valve

The actual position is shown with one or two digits on the LED display (1). For valve drives with two positions **L** for LOAD and **I** for die INJECT will be displayed.

Rear view of the Valve Drives S 6/-12/-16


Fig. 4 Rear view of the valve drives S 6/-12/-16

Except for the power socket 4 and an additional ground connector 3, the serial RS-232 interface 1 for external computer control and the remote control inputs and outputs 2 can be seen here.

Table 1 Terminal connections

GROUND:	Ground connector
HOME:	Sets the valve drive to position 1
FORWARD / LOAD:	Sets the valve drive to the LOAD position or the next higher position i.e. 2→3
BACKWARD / INJECT:	Sets the valve drive to the INJECT position or the next higher position i.e. 3→2
OUT:	Control output in position 1 has a low-resistance in the INJECT or in position 1 and has a high resistance in the LOAD or in positions 2-6 (2-16)

Depending on the DIP positions the four connectors **BIN 0** to **BIN 3** can be used alternatively for binary event control.

Operating the Valve Drive S 6/-12/-16

Power supply

The valve drive is equipped with a universal power supply unit which operates in a range of 90 to 260 Volts AC. A manual setting of the supply voltage is, therefore, not required.

Make sure that the main power supply is properly grounded and the correct power cable is used.

The power supply for the drive is not included in the housing. Please connect the power supply to the drive (item 4 in Abb. 4). Then connect the power cable to the power supply and the socket.

If required connect the ground connector 3 to the corresponding ground contact of any other **grounded** instrument.

Wait for at least 30 seconds before restarting the drive after switching it off. Otherwise an error message will be displayed.

Mounting the valve

The Smartline Electric Valve Drives S 6/-12/-16 have to be used in combination with one of the KNAUER-valves. For instructions on how to operate the valves refer to the corresponding valve manuals.

SOP 1 Mounting the valve to the electric valve drive

1. Fix the adapter plate with the two M4 screws at the rear side of the valve.
2. Set the drive on position 1 respectively LOAD.
3. Push the valve on the driving axle on the front side of the instrument.
4. Take care that the mark on the valve matches this one above the wheel. The valve position 1 must be on top.
5. Fix the valve and adapter plate with the two mounting screws.

If the position was not set on LOAD or 1, you will get the error message E 4 when starting the drive.

Configuring the Smartline Electric Valve Drive

At the bottom of the instrument the DIP-switches 1-4 (Abb. 5) are arranged. They are needed to adapt the drive to the actual mounted valve.


Fig. 5 DIP-switches 1-4 at the drive bottom

DIP 1 adapts the two position or multi-position valve (S 6).

The valve drive S 12, exclusively used to control the 13-port/1-channel switching valve, can only be operated with the DIP 1 switch in the ON position!

The valve drive S 16, exclusively used to control the 17-port/1-channel switching valve, can only be operated with the DIP 1 switch in the OFF position!

Using DIP 2 the control mode is set to either FORWARD/ BACKWARD/HOME or to the absolute position input. This will simultaneously activate the inputs BIN 0 through BIN 2.

The connector OUT/BIN3 can be used either as output or as input. The selection is performed by setting DIP 3 and DIP 4.

The DIP-switches 3 and 4 always should be set in opposite positions.

Not all of the possible combinations of DIP positions are permitted when operating the valve drive. If an improper DIP setting is used, an error message will appear either when switching the drive on or when the valve is switched.

The DIP setting always has to be done with the drive switched off. Otherwise the DIP positions will not become active.

Table 2 DIP-switch positions in valve drive S 6

DIP 1	ON	360° valve (7-p/1-ch)
	OFF	Two position valves (6-p/2-ch and 6-p/3-ch)
DIP 2	ON	binary control BIN 0-2 (360° valves only)
	OFF	RS-232 control/ manually/ HOME/ INJ./ LOAD
DIP 3 DIP 4	ON OFF	OUT/ BIN 3 output
DIP 3 DIP 4	OFF ON	not allowed

Table 3 DIP-switch positions in valve drive S 12

DIP 1	ON	360° valve (13-p/1-ch)
	OFF	not allowed
DIP 2	ON	binary control BIN 0-2
	OFF	RS-232 control/ manually/ HOME/ INJ./ LOAD
DIP 3 DIP 4	ON OFF	OUT/ BIN 3 output
DIP 3 DIP 4	OFF ON	OUT/ BIN 3 input

Table 4 DIP-switch positions in valve drive **S 16**

DIP 1	ON	not allowed
	OFF	360° valve (17-p/1-ch)
DIP 2	ON	binary control BIN 0-2
	OFF	RS-232 control/ manually/ HOME/ INJ./ LOAD
DIP 3 DIP 4	ON OFF	OUT/ BIN 3 output
DIP 3 DIP 4	OFF ON	OUT/ BIN 3 input

The data in the above tables imply that only one or two DIP combinations are possible for any of the valves to be controlled.




Table 5 Allowed DIP combinations

Valve	Mode	DIP position				Drive
6-port/2-channel 6-port/3-channel	analog	OFF	OFF	ON	OFF	S 6
7-port/1-channel	analog	ON	OFF	ON	OFF	
	binary	ON	ON	ON	OFF	
13-port/1-channel	analog	ON	OFF	ON	OFF	S 12
	binary	ON	ON	OFF	ON	
17-port/1-channel	analog	OFF	OFF	ON	OFF	S 16
	binary	OFF	ON	OFF	ON	
DIP switch No :		1	2	3	4	

6-Port-Multi Channel Valve Switch Positions :

The internal port connections of the multi-channel valves are shown in the following table.

Table 6 Port connections in multi-channel valves

Operated Key	 (INJECT)	 (LOAD)
Display	I	L
Connected ports 2-Channel valve		
Connected ports 3-Channel valve		

Controlling the Valve Drives S 6 /-12 /-16

Manual control

In order to operate the electrical valve drives S 6/-12/-16 manually, the DIP 2 must be in the OFF position.

2 position valves

The positions LOAD and INJECT can be set by directly using the  and  keys.

In the manual mode it is still possible to integrate the valve drive into a more complex chromatographic system. By switching the valve to the INJECT position, the OUT connector is short circuited. This signal can start time dependent gradient programs. The OUT connector is only active when the valve is in the INJECT position.

To avoid unintentionally restarting the program, turn the valve back to the LOAD position before the current program is completed.

6 , 12, and 16 position valves

Any valve position can be achieved in an increasing order by pushing the  key as often as necessary. The steps 6→1, 12→1, or 16→1 are also considered an increasing order. The opposite steps: 1→6, 1→12, or 1→16 are considered a decreasing order. To decrease a step press the  key.

Similar to the two position valves, the OUT connector can be used to control other instruments. The output is only active if the valve is in the HOME position (position 1). It is inactive if the valve is in positions 2 to 6, 2 to 12, or 2 to 16.

External Control of the Valve Drives S 6/-12/-16

The valve drive S 6/-12/-16 can be controlled with external contacts and can also send signals to other instruments. Communication is possible via the RS-232 interface or using the additional remote control connector (Pos 2 in Fig. 4).

If more than one valve drives are integrated into the system they should be controlled in the same control mode (RS-232 or WAGO). Specifically, if a coupled-columns technique is required.

Control via RS-232 interface

To control any valve via the RS-232 interface the drive has to be set to the analog mode (DIP switch 2 OFF).

The valve drives can be integrated into HPLC systems which are controlled with software packages such as ChromGate[®] or EuroChrom[®]. Instructions can be found in the corresponding software manuals.

The control via the serial RS-232 interface at the rear of the instrument, item 1 in Abb. 4 on page 20, also can be performed with small selfdesigned programs which can be integrated into non-Knauer software products.

The specifications for data transfer are:

9600 baud
8 bit
1 start-bit
1 stop-bit

The possible commands are listed in the following table. Any command must be confirmed using ↵ (= <ENTER>).

Table 7 Table of commands

Command	Action	Remark
1,2,3,...,6 (...,15,16,01,...)	Valve moves to the entered position	360°-Valves only
H,h	Valve moves to the HOME position (1)	360°-Valves only
L,l	Valve moves to the LOAD position	L/I -Valves only
I,i	Valve moves to the INJECT position	L/I -Valves only
V,v	Program version inquiry	
P,p	Actual position inquiry	
T,t	Actual valve type inquiry	

Short circuit signals from the OUT connector can be used besides the RS-232 control to control other instruments. The output is active or inactive under the same conditions as described for the manual control. If the additional inputs HOME, INJ., and/or LOAD are in use, the RS-232 connections have priority.

Remote control via Events

Generally the connections on the remote control strip are layed out as open collector (OC). Therefore the right polarity of any connection has to be obeyed.

2 position valves

Using this valve, the HOME input is inactive. If the LOAD input is short circuited to the GROUND, only the INJ. and GROUND have to be connected to the signal sending instrument (i.e. the interface box). If no signal is sent to the INJ., the valve will be in the LOAD position due to the LOAD GROUND short circuit. If INJ. is short circuited as well, the valve will switch to the INJECT position. INJ. has a higher priority than the LOAD.

6, 12, and 16 position valves

The 360° valves only can be operated either in an analog or in a binary mode. For the corresponding DIP positions refer to table 4 on page 10.

Controlling in analog mode

Short circuit signals from the OUT connector can be used to control other instruments. The output is only active if the valve is in the HOME position (1). It is inactive with the valve in positions 2 to 6, 2 to 12, and 2 to 16.

Three input connectors are available. A short circuit signal to HOME moves the valve to position 1. Similar signals to BACKW./INJ. and FORW./LOAD will move the valve back or forward one position. To switch the valve twice in the same direction the short circuit has to be shortly interrupted like a pulse circuit.

The most frequently used connector for the remote control strip is the Knauer interface box IF2. This wiring is illustrated in Fig. 6. If the drive control takes place via a pump, refer to the corresponding pump manual for the wiring.


Fig. 6 Wiring the valve drive to the interface box IF2

Controlling in the binary mode

Switching the control into the binary mode will activate the inputs BIN 0 through BIN 2. The output OUT is active in positions 1 through 6, however it is inactive in the process of switching. This activity change can be used as an outgoing signal to control other instruments.

To operate the 17-port/1-channel valve, BIN 3, also defined as an input, is connected so that positions 9 through 16 can also be binary controlled. When operating this valve, an output signal cannot be sent.

Every position now is binary defined as shown in the following table.

Table 8 Binary defined positions

Pos.	BIN 0 ($2^0=1$)	BIN 1 ($2^1=2$)	BIN 2 ($2^2=4$)	BIN 3 ($2^3=8$)	binary result
1	0	0	0	0	0
2	1	0	0	0	1
3	0	1	0	0	2
4	1	1	0	0	3
5	0	0	1	0	4
6	1	0	1	0	5
7	0	1	1	0	6
8	1	1	1	0	7
9	0	0	0	1	8
10	1	0	0	1	9
11	0	1	0	1	10
12	1	1	0	1	11
13	0	0	1	1	12
14	1	0	1	1	13
15	0	1	1	1	14
16	1	1	1	1	15

Assembling a remote control cable

If the remote control is performed not using the RS-232 interface a different type of control cable must be assembled for each application. Use the set plug strips with connectors enclosed in the delivery to prepare your cable according the following SOP.

SOP 2 Connecting the plug strips

Use this SOP for connecting wires to the plug strips.

1. Press down the cable clip as shown in **Fehler! Verweisquelle konnte nicht gefunden werden.** using a small screw driver or the WAGO liner tool.
2. Insert the uninsulated end of the cable into the opening under the latch.
3. Remove the tool from the plug. The cable is now firmly anchored in the plug strip.


Fig. 7 Connecting the Plug Strips

Valve drives S 6/12/16 as parts of the Smartline Assistant 6000

The valve drives S 6/12/16 as well as the UV-Detector 200 and/or the Pump S 100 can be an integrative part of the Smartline Assistant 6000. Up to three of these instruments can be placed in one Smartline housing. Thus, even a whole HPLC system can be put into one modul.

Fig. 8 Front view of the Smartline Assistant 6000 (A5003V065), with:
 1 UV Detector 200,
 2 Smartline Pump 100 and
 3 injection valve (manual)

An electrically driven injection valve can be placed in the rear position instead of the fractionation valve as in Fig. 9.

The preparative Modul 6100 is another example. It is equipped with a Smartline Pump 100 and zwei two valve drives. The pump and the front valve are forming the injection system. The second valve drive in the rear is for fractionation and solvent recycling.

Fig. 9 Frontview of the Prepmodule 6100

- 1 Smartline Pump 100
- 2 Injection valve and
- 3 Fractionation valve

The GPC Sample Preparation Unit 6500 is another example. It consists of a Smartline Pump 1000 and two versions of the Smartline Assistant 6000. Five valves and one UV Detector 200 are integrated into the Assistant versions V416 and V044.

Fig. 10 Front view of the GPC Cleanup Unit 6500

The control connectors of all three instruments are on the rear side of the modul. The DIP switches for configuring the valve drive are according to the single instruments at the bottom.

Troubleshooting

Error messages will be displayed in the form of „En“ (n = error code). A question mark “?” will appear if an invalid command is entered.

Table 9 List of error codes, causes, and solutions

code	probable cause	solution
E0	valve has not switched	exchange valve sealing
E1	motor power consumption to high	exchange valve sealing
E2	position changing to slow	exchange valve sealing
E3	DIP 3 and 4 in wrong setting	correct the settings
E4	position 1 no longer corresponds to HOME	exchange valve sealing
E5	DIP 1 and 2 in wrong setting	correct the settings
	wrong binary code	correct the settings
E6	error while reading the EEPROM	restart the instrument

In case the purposed solutions will not solve your problems, please contact our service department.

Spare parts and accessories

Power supply cable, 230V	M1642
RS-232 cable	A0895
Set WAGO connecting rails	A1402
Set signal conductors	A1467

Maintenance by the user

The Electric Valve Drives S 6, S 12, and S 16 do not need any maintenance. In case of problems please contact our service department.

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Guarantee statement

The guarantee period of the Smartline Electric Valve Drives S 6, S 12, and S 16 is 12 months beginning from the date of dispatch from Berlin. Operation inconsistent with manufacturer's instructions or damage caused by unauthorized service personnel are excluded from guarantee. Damage caused by blockages and wear and tear parts such as fuses and seals are not covered by the guarantee. Claims under this guarantee are valid only if the enclosed guarantee card is returned to us at the address shown below within 14 days of receipt of the instrument. Defective valve drives should be sent to the manufacturer for repair.

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If we find a defect covered by the guarantee, repair or replacement, at our discretion, will be carried out free of charge. Packing and transport costs are borne by the purchaser.

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Declaration of conformity

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Manufacturer's name and address:
Herstellername und -adresse

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Smartline Valve Drives S6, S12, S16
with different valves

Order Numbers, Bestellnummern: **A 550xx, A551xx** and **A 552xx**

complies with the following requirements and product specifications:

- Low Voltage Ordinance (73/23/EWG)
EN 61010-1 (08/2002)
- Engineering Guidelines (89/392/EWG)
- EMC Ordinance (89/336/EWG)
EN 6100-3-2 (10/2006)
EN 61326-1 (10/2006)

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Dr. Ing. Herbert KNAUER GmbH