



BlueShadow[®] Assistant

Erläuterungen zu den Aufklebern auf der Geräterückseite

Auf der Geräterückseite sind Aufkleber mit Informationen über die Module und den Assistenten angebracht. Unterhalb des Barcodes stehen folgende Informationen:

- Seriennummer
- MAC-Adresse

Abbildung 1 zeigt die Position der Aufkleber auf der Rückseite. Abbildung 2 zeigt die Position der Module in der Aufsicht.

Beschreibung	Abbildung
<ul style="list-style-type: none"> ① Aufkleber des Assistenten ② Aufkleber des Gerätes an Position A ③ Aufkleber des Gerätes an Position C ④ Aufkleber des Gerätes an Position B ⑤ Gerätevorderseite ⑥ Geräterückseite 	 <p>Abb. 1 Geräterückseite</p>  <p>Abb. 2 Anordnung der Module</p>



BlueShadow[®] Assistant

Comments on the Labels on the Rear Panel

The labels on the rear panel contain information about the modules and the assistant. The information printed below the bar code include the following:

- Serial number
- MAC address

Figure 1 shows the positions of the labels on the rear panel. Figure 2 shows the positions of the modules in plan view.

Description	Figure
<p>① Label of assistant</p> <p>② Label of module in position A</p> <p>③ Label of module in position C</p> <p>④ Label of module in position B</p> <p>⑤ Front panel</p> <p>⑥ Rear panel</p>	 <p>Abb. 3 Rear panel of the device</p>  <p>Abb. 4 Positions of the modules</p>