









Overview of 1/8" fittings for V 4.1 valves





The ports of V 4.1 valves for 1/8" capillaries are coned with a UNF 1/4-28 thread. Appropriate fittings are provided with each valve. The tables below give an overview on the available fittings. You can also find this overview on our website: www.knauer.net.

Bushings for 1/8" UNF 1/4-28 coned








Description	Image	Material	Installation type	max. back-pressure [bar]	Required sealing	Article number
Bushing, short		SST	wrench tight 5/16" Hex	300	Ferrule A7213	A7205
Bushing, long		SST	wrench tight 5/16" Hex	300	Ferrule A7213	A7206
Bushing, long		SST	wrench tight 5/16" Hex	300	Biconical seal ring A7216 or clamp ring A7214 + seal ring A7217	A7207
Blind fitting		SST	wrench tight 5/16" Hex	300	None	A7208
Bushing with integrated sealing ring (one-piece)		PEEK	Fingertight	50	None	A7209
Bushing		PEEK	Fingertight	100	Biconical seal ring A7216 or clamp ring A7214 + seal ring A7217	A7210
Bushing with integrated sealing ring (one-piece)		PCTFE	Fingertight	69	None	A7211

Description	Image	Material	Installation type	max. back-pressure [bar]	Required sealing	Article number
Blind fitting		PEEK	wrench tight 5/16" Hex	100	None	A7212

Ferrules, seal rings and clamp rings for 1/8" UNF 1/4-28 coned

Description	Image	Material	Article number
Ferrule		SST	A7213
Clamp ring		SST	A7214
Clamp ring		Titan	A7215
Biconical seal ring		PETP	A7216
Seal ring		PEEK	A7217

Adapters and couplings for 1/8" UNF 1/4-28 coned

Description	Image	Material	Connection 1	Connection 2	Article number
Coupling for 1/8" <-> 1/8" capillary		SST	Internal thread 1/8" UNF 1/4-28 coned	Internal thread 1/8" UNF 1/4-28 coned	A7218
Coupling for 1/8" <-> 1/8" capillary		Titan	Internal thread 1/8" UNF 1/4-28 coned	Internal thread 1/8" UNF 1/4-28 coned	A7219
Coupling for 1/8" <-> 1/8" capillary		PEEK	Internal thread 1/8" UNF 1/4-28 coned	Internal thread 1/8" UNF 1/4-28 coned	A7220
Coupling for 1/16" <-> 1/8" capillary		SST	Internal thread 1/16" UNF 10-32 coned	Internal thread 1/8" UNF 1/4-28 coned	A7221
Coupling for 1/16" <-> 1/8" capillary		Titan	Internal thread 1/16" UNF 10-32 coned	Internal thread 1/8" UNF 1/4-28 coned	A7222
Coupling for 1/16" <-> 1/8" capillary		PEEK	Internal thread 1/16" UNF 10-32 coned	Internal thread 1/8" UNF 1/4-28 coned	A7223
Adapter 1/8" to 1/16"		SST	Internal thread 1/16" UNF 10-32 coned	External thread 1/8" UNF 1/4-28 coned	A7204
Adapter 1/8" to 1/16"		PEEK	Internal thread 1/16" UNF 10-32 coned	External thread 1/8" UNF 1/4-28 coned	A7224

Analytical
HPLC

Multi-Column
Chromatography,
SMB

Preparative
HPLC

FPLC

Osmometry

Dosing,
Metering,
Pumping

Detection

Technical data are subject to change without notice.


All trademarks are the property of their respective owners.
Our general terms and conditions apply:
www.knauer.net/terms

© KNAUER 2020 V6886
Printed on environmentally friendly paper.

KNAUER

Wissenschaftliche Geräte GmbH
Hegauer Weg 38
14163 Berlin, Germany

+49 30 809727-0
+49 30 8015010 (Fax)
info@knauer.net
www.knauer.net



KNAUER